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CMV RQ Kit Manual

**For Real-Time PCR Detection and
Quantitation of Cytomegalovirus
(CMV) DNA**

For use with Rotor-Gene or StepOne
Research use only

NG-WI-ASL-05-400
Version 4.0
March 2017


Table of Contents:

1. Introduction.....	2
2. Kit Contents	3
3. Storage and Stability.....	3
4. General Precautions	4
5. Additionally Required Materials.....	4
6. Specimen, storage and transport	5
7. Interfering substances.....	5
8. Internal control (IC).....	5
9. DNA isolation	5
10. PCR Protocol	6
11. Programming of the Rotor-Gene	6
12. Programming of StepOne(/Plus)	7
13. Programming other machines	7
14. Data Analysis: Rotor-Gene.....	8
15. Data Analysis: StepOne	10
16. Quantitation	12
17. Linear Range	12
18. Sensitivity	12

CMV RQ kit is intended for use with Rotor-Gene or StepOne machines and for the quantitative detection of CMV-DNA. This kit is for research use only.

1. Introduction

Human Herpesvirus 5 (HHV-5) known as Cytomegalovirus (CMV) is a double-stranded DNA virus and a member of Herpesviridae family. It is a widespread pathogen with high seropositivity (about 80%) in adult population. Primary infection is usually subclinical and mostly happens in early childhood. The primary infection is followed by the latent infection and CMV persists in the infected person through life. CMV can be reactivated in immunocompromised hosts for example after stem cell, bone marrow or solid organ transplantation or in the course of HIV infection. In these patients, viral load assays are the cornerstone for the diagnosis and monitoring of CMV disease. Viral load quantification in blood allows assessing patients at risk of developing relapsing infection, early detection of active CMV infection, guiding preemptive therapy, monitoring response to therapy, determining duration of therapy and prediction of treatment failure due to emergence of resistant strains.

CMV RQ kit provides a ready-to-use Real-Time PCR system for detection and quantitation of cytomegalovirus with Rotor-Gene or StepOne machines. This method provides the highest sensitivity and widest dynamic range among other available methods. In this method application of fluorescent dye labeled probes allows detection of amplified product.

Analysis of fluorescent kinetics also leads to quantification of the target sequence in the reaction without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product.

This kit also incorporates an *Internal Control* (IC) to identify possible PCR inhibition.

2. Kit Contents

The kit contains a manual, a CD with Rotor-Gene and StepOne templates and following reagents:

Label	Content	Quantity
CMV Mix	PCR Master mix*	360 μ l
CMV S1	Standard 1: 10,000 copy/ μ l	150 μ l
CMV S2	Standard 2: 1,000 copy/ μ l	150 μ l
CMV S3	Standard 3: 100 copy/ μ l	150 μ l
CMV S4	Standard 4: 10 copy/ μ l	150 μ l
Water	PCR Grade Water	200 μ l

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits, respectively.

3. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws more than few times to prevent reduced sensitivity.

4. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components **on ice** completely, mix by flickering followed by a quick spin and **store on crushed ice** after.
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

5. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Table top microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.7ml microtubes and PCR microtubes
- Disposable powder-free gloves

6. Specimen, Storage and Transport

Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or citrate plasma for CMV detection. Whole blood or plasma should be shipped at +4°C. Upon receipt plasma should be separated from whole blood and can be stored at +4°C for few days or aliquoted and stored at -20°C for up to few weeks.

7. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes should not be used. Samples of heparinized patients must not be used as well. Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

8. Internal Control (IC)

To examine the possible PCR inhibition and to prevent false negative results, the kit contains an *Internal Control* included in PCR Mix. Internal control should generate a CT of 28-32 in Yellow Channel on Rotor-Gene and a CT of 28-34 in VIC channel on StepOne.

9. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)

- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)
- QIAampUltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAampMiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

10. PCR Protocol

Thaw the reagents on ice completely followed by a brief mixing and a quick spin. Place required number of tubes on cold block. Consider one tube for each sample plus one for each standard and one for the negative control.

Pipette 15ul of CMV Mix directly to each tube followed by adding 10ul of standard or isolated DNA.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring.

11. Programming Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the CD provided in the kit and double click on CMV Template 0.2 or CMV 0.1 according to the tubes used. Program starts. Click on Start button (Green button on the top menu). On

the pop up window click start again and save program on desired location.

12. Programming StepOne

Open the StepOne software (V 2.*). On the Set Up menu click on Template and select the file on CD provided with the kit. Click on Plate Setup. One negative control, 4 standards and 10 samples are defined. You may change plate set up using right click options (copy, paste, clear). You may also add/remove samples or change sample names on "Define Targets and Samples" menu. When finished, click on "Start Run" and save the experiment on desired location. Instrument will start shortly.

13. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 10 min	1
2	95°C x 15 sec	45
	58°C x 60 sec	

Fluorescence should be collected at 58°C for FAM and VIC dyes. The CMV Mix contains ROX. Final concentration of ROX in reaction is 300nM.

14. Data Analysis: Rotor-Gene

Before analyzing results, make sure in the sample menu all the standards have been defined as “standard” and the relative concentrations have been entered. Patient samples should be defined as “unknown” and Negative control or no template control as “Negative Control” or “NTC” respectively.

Analyze data according to manufacturer recommendations. Perform quantitative analysis for **CMV (Green channel)** and qualitative analysis for **Internal Control (Yellow channel)**. Briefly, click on analysis menu and then under Quantitation tab double click on cycling A. Green.

In the pop up for Automatic Threshold increase the minimum or lower bound until it surpasses the negative control or the NTC fluorescence, and then click on OK.

Repeat the above for Cycling A. Yellow but cancel the Automatic Threshold and manually put threshold on 0.1. Figures 1 and 2 represent typical graphs for Rotor-Gene machine.


Figure 1. Typical CMV Graph in Green Channel for Rotor-Gene


Figure 2. Typical IC Graph in Yellow Channel for Rotor-Gene

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used.

In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.

Consider following points when analyzing:

- A sample is **Positive** if it is positive in Green channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results in the Cycling A. Green are valid.
- A sample is **Negative** if it is negative in Green channel while it is positive in Yellow channel with a sigmoid graph and CT of 28-32.

- Results are **Inconclusive** and the test should be repeated if a sample is negative in both of Green and Yellow channels.

15. Data Analysis: StepOne

Analyze data according to manufacturer recommendations. Briefly, click on Analyze and set the threshold for CMV/FAM at 0.1 and at 0.05 for IC/VIC.

Figures 3 and 4 represent typical graphs for StepOne machine.


Figure 3. Typical CMV Graph in FAM Channel for StepOne


Figure 4. Typical IC Graph in VIC Channel for StepOne

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used.

In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.

Consider following points when analyzing:

- A sample is **Positive** if it is positive in FAM/CMV channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results are valid.
- A sample is **Negative** if it is negative in FAM/CMV channel while it is positive in VIC/IC channel with a sigmoid graph and CT of 28-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both of FAM/CMV and VIC/IC channels.

16. Quantitation

The kit provides 4 quantitation standards with defined titers to generate a standard curve for quantification of samples viral load. Working with Rotor-Gene machine, the standard curve from a previous run can also be imported for quantification of samples to the recent run. To do so, at least one standard must be used in the current run. Apparently using all five standards in each run will lead to more accurate results.

Quantitation standards are defined as copy/ μ l. To convert the result to copy/ml following equation should be used:

$$\text{Result (copy/ml)} = \frac{\text{Result(copy/}\mu\text{l)} \times \text{elution volume (}\mu\text{l)}}{\text{sample volume (ml)}}$$

“Sample volume” is the plasma volume used for DNA isolation and “Elution volume” is the volume of buffer or water used to elute or dissolve isolated DNA.

17. Linear Range

The linear range of the kit was assessed with dilution series of the cloned target and showed to be linear in the range of 1,000,000 copies/ μ l to 1 copy/ μ l.

18. Sensitivity

The analytical detection limit of the kit was assessed with dilution series of the cloned target and showed a limit of detection equal to 1 copy/ μ l.