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TB RQ Kit Manual

For Real-Time PCR Detection and Quantitation of Mycobacterium Tuberculosis Complex

For use with Rotor-Gene or StepOne
Research use only

NG-WI-ASL-10-303

Version 3.3

June 2018

Table of Contents:

1. Introduction	2
2. Kit Contents	3
3. Storage and Stability	3
4. General Precautions	3
5. Additionally Required Materials.....	4
6. Specimen, storage and transport	4
7. DNA isolation	4
8. Internal control (IC).....	5
9. Quantitation	5
10. PCR Protocol	6
11. Programming of Rotor-Gene	7
12. Programming of StepOne(/Plus)	7
13. Programming of other machines	6
14. Data Analysis: Rotor-Gene.....	9
15. Data Analysis: StepOne	11
16. Linear Range	11
17. Sensitivity	12

TB RQ kit is intended for the quantitative detection of *Mycobacterium tuberculosis* complex DNA with Rotor-Gene or StepOne machines. This kit is for research use only.

1. Introduction

Tuberculosis (TB) is one of the major public health problems worldwide and results in about 2 million deaths annually. This infection is caused by a group of closely related species and subspecies of intracellular pathogens called *Mycobacterium tuberculosis* complex (MTC) including *M. tuberculosis*, *M. bovis* and *M. bovis* BCG.

While bacterial culturing is still the gold standard for TB diagnosis, Real-Time PCR is increasingly used as an alternative due to its rapid, sensitive and specific results. It also provides the highest sensitivity and widest dynamic range among other methods. In Real-Time PCR, application of fluorescent dye labeled probes allows detection of amplified product. Analysis of fluorescent kinetics also leads to quantification of the target sequence in the reaction without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product.

TB RQ kit provides a ready-to-use Real-Time PCR system for detection and quantitation of *Mycobacterium tuberculosis* complex DNA with Rotor-Gene (Qiagen) or StepOne (Applied Biosystems) machines. The kit also incorporates an *Internal Control* (IC) to identify possible DNA extraction failure or PCR inhibition or set up errors.

2. Kit Contents

The kit contains a manual, a CD with Rotor-Gene and StepOne templates and following reagents:

Label	Content	Quantity
TB Mix	PCR Master mix*	360 µl
TB S1	Standard 1: 100,000 copies/µl	150 µl
TB S2	Standard 2: 10,000 copies/µl	150 µl
TB S3	Standard 3: 1,000 copies/µl	150 µl
TB S4	Standard 4: 100 copies/µl	150 µl
TB S5	Standard 5: 10 copies/µl	150 µl
TB IC	Internal Control*	250 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

3. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws more than few times to prevent reduced sensitivity.

4. General Precautions

In order to prevent false results, always pay attention to the following points:

- Treat all samples as potentially infectious.
- Decontaminate samples prior to DNA extraction.
- Within the pre-PCR work area assigns three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the TB Mix is aliquoted into tubes; and c) Reaction preparation area for addition of extracted DNA to the tubes.

- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working**.
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

5. **Additionally Required Materials**

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Table top microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipetters and nuclease free filtered tips
- Mycobacterium/TB DNA extraction kit
- Nuclease free 1.7ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

6. **Proper Specimen**

Sputum, BAL (bronchoalveolar lavage), bronchial secretion, CSF and peritoneal puncture are proper patient samples.

7. **DNA Isolation**

DNA isolation can be performed using different kits from various manufacturers. Make sure the kit is especially manufactured for Mycobacterium/TB DNA extraction. If not available, you may use QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany) with some modification. For more details

please refer to "artus® M. tuberculosis RG PCR Kit Handbook (Mat# 1046960EN, Qiagen)".

To monitor DNA extraction process, internal control should be applied to the extraction process. For more details, please, refer to section 8 of this handbook.

8. Internal control (IC)

In order to evaluate the possibility of DNA extraction failure and PCR inhibition and prevent false negative results, TB RQ kit contains internal control (IC). IC can be used during extraction process or simply added to TB Mix.

To monitor of both DNA extraction and PCR reaction, internal control should be added to the mixture of lysis buffer and patient sample. Required volume of internal control is 10% of elution buffer. For instance if extracted DNA is eluted with 100ul, then 10ul of internal control should be added to the mixture of lysis buffer and patient sample. **Please note that internal control should not be added directly to the patient sample (i.e. before addition of lysis buffer) as it loses its efficiency.**

If IC is added to TB Mix, only PCR inhibition can be monitored. For this purpose, 1 microliter of IC reagent should be added to each reaction. For example, for 10 PCR reactions, 15ul of internal control should be added to 150ul of TB Mix before added to the tubes.

In a successful DNA extraction and PCR test, IC should generate a CT of 31-37 in Yellow/VIC Channel.

9. PCR Protocol

Thaw the reagents on ice completely followed by a brief mixing and a quick spin. Place required number of tubes on cold block.

Consider one tube for each sample plus one for each standard and one for the negative control.

If IC is introduced during extraction process, pipette 15µl of TB Mix directly to each PCR tube.

If IC is added to the TB Mix, add 15ul of the [prepared mix](#) (as described in section 8) to each PCR tube.

Then add 10ul of extracted DNA, [standard](#) or water to each tube.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring.

10. Programming Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the CD provided in the kit and double click on TB 0.2 V.3 if you are using 0.2ml tubes. If you are using 0.1ml strip tubes, double click on TB 0.1 V.3. Program starts. Click on Start button (Green button on the top menu). On the pop up window, click start again and save program on desired location.

11. Programming StepOne

Open the StepOne software (V 2.*). On the Set Up menu click on Template and select the file on CD provided with the kit. Click on Plate Setup. One negative control, 5 standards and 10 samples are defined. You may change plate set up using right click options (copy, paste, clear). You may also add/remove samples or change sample names on "Define Targets and Samples" menu.

When finished, click on “Start Run” and save the experiment on desired location. Instrument will start shortly.

12. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 10 min	1 cycle
2	95°C x 15 sec	40 cycles
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. The TB Mix contains ROX. The final reaction concentration of ROX is 300nM.

13. Data Analysis: Rotor-Gene

Before analyzing results, make sure in the sample menu all the standards have been defined as “standard” and the relative concentrations have been entered. Patient samples should be defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC", respectively.

Analyze data according to manufacturer recommendations. Perform quantitative analysis for **TB (Green channel)** and qualitative analysis for **Internal Control (Yellow channel)**. Briefly, click on analysis menu and then, under Quantitation tab double click on “Cycling A. Green”.

In the pop up for Automatic Threshold increase the minimum or lower bound until it surpasses the negative control or the NTC fluorescence, and then, click on OK.

Repeat the above for “Cycling A. Yellow” but, cancel the Automatic Threshold and manually put threshold on 0.2.
 Figures 1 and 2 represent typical graphs for Rotor-Gene machine.


Figure 1. Typical TB Graph in Green Channel for RotorGene


Figure 2. Typical IC Graph in Yellow Channel for RotorGene

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used.

In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.

Consider following points when analyzing:

- A sample is **Positive** if it is positive in Green channel with a sigmoid graph and a CT of less than 35. The titer and quantitation results in the Cycling A. Green are valid.
- A sample is **Negative** if it is negative in Green channel while it is positive in Yellow channel with a sigmoid graph and CT of 28-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both of Green and Yellow channels.

14. Data Analysis: StepOne

Analyze obtained data according to manufacturer recommendations. Briefly, click on “Analyze” and set the threshold for **TB/FAM** at 0.1 and at 0.1 for **IC/VIC**.

Figures 3 and 4 represent typical graphs for StepOne machine.


Figure 3. Typical TB Graph in FAM Channel for StepOne


Figure 4. Typical IC Graph in VIC Channel for StepOne

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used.

In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.

Consider following points when analyzing:

- A sample is **Positive** if it is positive in FAM/TB channel with a sigmoid graph and a CT of less than 35. The titer and quantitation results are valid.
- A sample is **Negative** if it is negative in FAM/TB channel while it is positive in VIC/IC channel with a sigmoid graph and CT of 28-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both of FAM/TB and VIC/IC channels.

15. Quantitation

The kit provides 5 quantitation standards with defined titers to generate a standard curve for quantification of samples mycobacterial load. Working with Rotor-Gene machine, the standard curve from a previous run can also be imported for quantification of samples to the recent run. To do so, at least one standard must be used in the current run. Apparently using all five standards in each run will lead to more accurate results.

Quantitation standards are defined as copy/ul. To convert the result to copy/ml following equation should be used:

$$\text{Result (copy/ml)} = \frac{\text{Result (copy/}\mu\text{l)} \times \text{elution volume (}\mu\text{l)}}{\text{sample volume (ml)}}$$

“Sample volume” is the sample volume used for DNA isolation and “Elution volume” is the volume of buffer or water used to elute or dissolve isolated DNA.

16. Linear Range

The linear range of the kit was assessed with dilution series of the cloned target and showed to be linear in the range of 1,000,000 copies/ul to 10 copies/ul.

17. Sensitivity

The analytical detection limit of the kit was assessed with dilution series of the cloned target and showed a limit of detection equal to 1 copy/ul.