

LTEK INNO & INNO-M

Microplate Spectrophotometer

Absorbance & Luminescence

New generation for Microplate spectrophotometers is here



Microplate Spectrophotometers

INNO (Absorbance) & INNO-M (Absorbance + Luminescence)

The LTEK™ INNO microplate spectrophotometer is a 21st century new generation instrument for all the researches that deals with the absorbance and luminescence. INNO and INNO-M are based with monochromatic measuring technology which allows your research to higher level. INNO-X which is the software for INNO and INNO-M has powerful functions with high accuracy of result data.

One of the most convenient part of INNO-X is that the end users can download

the software from the website of LTEK at anytime and anywhere only with the internet connection. And be able to update the software whenever the new upgraded version is available.

Allowing variety of endpoint, kinetic, spectral scanning, and area scanning assays. Results are provided quickly and presented conveniently on Microsoft excel which allows you to be able to save and manage the result data securely.

LTEK INNO Microplate Spectrophotometers come in two different versions INNO & INNO-M.

INNO Microplate

Spectrophotometer Absorbance.

- 200 ~ 1000nm detectable wavelength range with 1nm increment.
- · Wavelength selection monochromator.
- · Xenon flash lamp for semi-permanent life time.



INNO-M Microplate

Spectrophotometer Absorbance & Luminescence.

- \cdot 200 ~ 1000nm detectable wavelength range with 1nm increment.
- · Wavelength selection monochromator.
- · Xenon flash lamp for semi-permanent life time.



INNO-X and the performance

Able to read 200 ~ 1000nm without filters

Unlike other basic reader instruments out there, INNO & INNO-M use monochromator system instead of filters with a proven optical design that helps ensure excellent sensitivity and unparalleled results. Brining wide range of wavelengths from 200 to 1000nm with just 1 nm increment.

Quantitative and Qualitative measuring protocol option screen



7 different regressions to view the result data in Assay and Calibration.

Quantitative and Qualitative assay kit test

Various types of Quantitative and Qualitative assay tests are available. In each protocol there are 7 different regressions that can be used for the analyzing data.


Linear, Quadratic, Cubic, Log, Exponential, point-to-point, 4PL

Freely utilizable data on PC

Save the protocol for convenient use every time.

Access to the well organized and high-quality

User interface to obtain new generation result data.







Multiple options to maximize measurement performance.


With End Point / Kinetic settings, up to 6 wavelengths are possible to be set.

In Spectral scanning from 200 to 1000nm can be read with 1nm increment.

Area scanning function read up to 338 points per well.

6 ~ 384well plate read

Variety types of plates (6 \sim 384 well plates) can be read.




24well DNA/RNA Quantitative measurement


Using 2ul DNA / RNA samples Quantitative values are measurable.


Unknown quantitative values of DNA / RNA samples can be read from $240 \sim 320$ nm to check the quantitative values.

With 24 wells supporting function, variety samples (DsDNA, RNA, ssDNA, 1Abs at 1cm = 1mg/ml BSA, IgG and Lysozyme) concentration value can be measured.















- Monochromator-based optical system for free selection of wavelengths from 200nm to 1000nm.
- · No filter is needed for these readings.
- · INNO is able to read microplate with 6 ~ 384 wells
- INNO can perform endpoint, Kinetics reading, and spectral scanning• Both photometric acuity and linearity of INNO-M, it should be 0-2,000 OD +/- 1%.
- INNO is able to be used in all studies such as routine biology tests, protein analysis, nucleic acidity.
- Nano-V can analyze DNA/RNA quantities. (Micro volume plate / 24well supported)
- · Plate capable of quantitative analysis with 2ul.
- · The light source of INNO is the Xenon lamp.
- The software to be supplied with INNO-M; Abs, UV-Abs, and modes, Endpoint, Kinetic and Spectral scanning.
- · INNO-M is able to report result graphics in excel file format.
- INNO-M software supports Linear, Quadratic regression, Cubic regression, Log regression, Exponential regression, Linear logarithmic regression, point-to-point, and 4PL regression graphic options.





· Mea

Luminescence

- · Detection method is Photomultiplier (PMT)
- · Measuring Wavelength range is between 300nm to 700nm.
- · Limit of Detection 3x10⁻²¹ moles.



LTEK

Specification for INNO and INNO-M

Absorbance	
Wavelength Accuracy	±2 nm
Electrical Requirements	INPUT 100 to 240V 50/60Hz / (65W Adaptor)
Microplate type	6 ~ 384 well plate
Detector	Photodiode
Light source	Xenon flash
Wavelength Range	200 to 999 nm
Wavelength selection	Monochromator
Application	Wavelength scanning, end point, Kinetic, Area scan
Dynamic range	0 ~ 4.0 OD
OD accuracy	0 ~ 2 OD ±1%
OD linearity	0 ~ 2 OD ±1%
OD repeatability	0 ~ 2 OD ±1%
Shaking	Two step speed
Software	INNO X (Windows Software)
DNA/RNA Micro Volume plate	24well / 2ul Sample (Option)
LUMINESCENCE	
Detector	Photomultiplier (PMT)
Wavelength range	300 – 700nm
Peak Wavelength	420nm
Limit of Detection (moles)	3x10 ⁻²¹ moles
Supported software regression	Linear, Quadratic, Cubic, Log, Exponential, Linear logarithmic, point- to-point, 4PL
Weight	8kg
Size	333x303x245

Available products

Product name	Description
INNO	Microplate Spectrophotometer (Absorbance)
INNO-M	Microplate Spectrophotometer (Absorbance + Luminescence)
INNO-I	1.5ml Microcentrifuge Cuvette Type Luminometer (Luminescence)
INNO-N	Micro Volume Spectrophotometer
NANO-V(Option)	24 Wells / 2ul samples

Contact info

Homepage: www.ltekc.com / Tell: +82-70-7755-9375 / Overseas sales manager: ysk0109@ltekc.com

