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
# **Factor II RG Kit Manual**

**For Real-Time PCR Detection of  
FactorII (G20210A) mutation**

For use with Rotor-Gene 6000  
(Corbett Research)

Research use only  
NG-WI-ASL-15-101  
Version 1.2  
Desamber 2019

## Table of Contents:

1. Introduction.....	2
2. Kit Contents .....	3
3. Storage and Stability .....	3
4. Additionally Required Materials .....	3
5. General Precautions .....	4
6. Specimen, storage and transport .....	4
7. Interfering substances .....	5
8. DNA isolation .....	5
9. Internal Control (IC) .....	5
10. PCR Protocol .....	5
11. Programming of Rotor-Gene6000 .....	6
12. Data Analysis.....	6

**Factor II RG kit** is intended for the detection of prothrombin G20210A mutation in human DNA with Rotor-Gene machine. This kit is intended for research use only

## 1. Introduction

Prothrombin or factor II is the precursor for thrombin in the coagulation cascade. G20210A mutation in prothrombin gene is the second most prevalent genetic risk factor for venous thromboembolism (VTE) after Factor V Leiden. This mutation increases the risk of venous thrombosis about 2 to 7 times. It has also been linked to increased risk of pregnancy loss. Prevalence of this mutation ranges between 0.5 to 5% depending on the race and ethnic background.

The diagnosis of factor II G20210A mutation requires DNA analysis of *F2*, the gene encoding factor II.

F2 RG kit provides a ready-to-use Real-Time PCR system for detection of G20210A mutation with Rotor-Gene machine. This method applies fluorescent probes allowing detection and analysis of amplified product while the reaction is in progress. Provided fluorescent kinetics also leads to assessment of the target sequence and no further post-amplification analysis is required. Therefore, turn-around time is reduced and the possibility of contamination with the PCR product is removed.

## 2. Kit Contents

The kit contains a manual, a CD with Rotor-Gene template and following reagents:

Label	Content	Quantity
F2 Mix	PCR Master mix*	480 $\mu$ l
MM Pos Ctrl	MM Positive Control	50 $\mu$ l
MW Pos Ctrl	MW Positive Control	50 $\mu$ l
WW Neg Ctrl	WW Negative Control	50 $\mu$ l
Water	PCR Grade Water	200 $\mu$ l

\*1, 2 or 4 tubes for 24, 48 or 96 reaction kits.

## 3. Storage and Stability

All the kit components should be shipped and stored at  $-20^{\circ}\text{C}$  and are stable until the expiry date mentioned. Avoid repeated freeze-thaws more than 3 times to prevent reduced sensitivity.

## 4. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Table top microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tip
- DNA extraction kit
- Nuclease free 1.7ml microtubes and PCR microtubes
- Cold block

## 5. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area assign three separate spaces for a) Sample storage and extraction, b) Reagent preparation where the mastermix is aliquoted into tubes and c) Reaction preparation area for addition of extracted sample to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- **Thaw kit components on crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice after.**
- Do not place 0.2ml PCR tubes on crushed ice. Use cold blocks instead.

## 6. Specimen, Storage and Transport

Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or citrate as anticoagulant for Factor II mutation detection. Whole blood should be shipped and stored at +4°C (stable for 3 days). Whole blood or buffy coat can be stored at +4°C for few days. Otherwise should be aliquoted and stored at -20°C which is stable for few months.

DNA can directly be extracted from whole blood or buffy coat.

## 7. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes should not be used. Samples of heparinized patients must not be used as well. Elevated levels of bilirubin ( $\leq 4.5$  mg/dl) and lipids ( $\leq 1000$  mg/dl) and hemolytic samples do not influence the extraction and PCR.

## 8. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend using following:

- High Pure PCR Template Preparation Kit (Cat. no. 11796828001, Roche Applied Science, Mannheim, Germany).
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany).

## 9. Internal Control (IC)

Each person carries wild type, mutant or both alleles of F2. These alleles serve as both the target and the Internal Control of the assay. Any sample should always be positive for at least one of them. If a sample is negative for both alleles, the reaction is inhibited and the test should be repeated.

## 10. PCR Protocol

Thaw the reagents on crushed ice completely followed by brief mixing and a quick spin. Place required number of tubes on

cold block. Consider one tube for each sample plus one for each control and one for NTC.

**Pipette 20µl of F2 Mix directly to each tube followed by addition of 5µl of isolated DNA.**

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring too.

## 11. Programming Rotor-Gene

*- Before you start the machine, make sure you have attached the ring on the rotor!*

*- Please note that at temperatures of 30°C and above, Rotor-Gene cooling system fails and test results are not valid! The optimum temperature for proper function of this machine is 20-25°C!*

Open the CD provided in the kit and double click on Factor II 0.2 RG or Factor II 0.1 RG (according to the used tubes) to open the program. Click on Start button (Green button on the top menu). On the pop up window click start again and save program.

On the pop up sample window, name the samples according to the list.

You may also set Rotor-Gene according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1

2	<b>95°C x 15 sec</b>	40
	<b>55°C x 30 sec*</b>	
	<b>72°C x 15 sec</b>	
3	<b>95°C x 15 sec</b>	1
4	<b>40°C x 10 min</b>	1
5	<b>Melt 42°C-75°C *</b>	1

Fluorescence should be collected at steps marked with star (\*) for FAM.

The PCR Mix is ROX free.

## 12. Data Analysis

Analyze data according to manufacturer recommendations.

Briefly, on analysis menu click on Melt tab and then double click on Melt A Green.

On the right side menu, uncheck the box for “Flip sign of dF/dT”. On the bottom of “Melt Curve Analysis” window click on “Auto-Scale” tab. Set threshold manually above the background, so only two picks remain above the threshold. **The left pick (at 48-49°C) is for wild type allele (W) and the right one (at 58-59°C) is for the mutant allele (M).** Analyze the results as follows:

- A sample with only a single pick at 48-49°C is **Negative** for Factor II G20210A and is considered **wild type homozygous (WW)**.
- A sample with only a single pick at 58-59°C is **Positive** and **homozygous (MM)** for Factor II G20210A (both alleles are mutant).
- A sample with two picks at 48-49°C and at 58-59°C is **Positive** for Factor II G20210A and **heterozygous (WM)**. Such a patient carries both wild type and mutant alleles for Factor II.

Results are inconclusive and the test should be repeated if:

- No pick is detected for a patient at either 48-49°C or 58-59°C.
- A pick is detected for NTC.

Genotypes can be called automatically if bins are defined for W and M picks. Please refer to the Rotor-Gene manual for further information. Figures 1, 2 and 3 represent sample analysis and the relevant peak for each genotype.


Figure 1: Wild type (WW) graph on Rotor-Gene



Figure 2: Mutant Heterozygote (WM) graph on Rotor-Gene

**Figure 3:** Mutant Homozygote (MM) graph on Rotor-Gene