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JC RQ Kit Manual

**For Real-Time PCR Detection and Quantitation of
JC Virus**

For use with Rotor-Gene or StepOne
Research use only

NG-WI-ASL-26-100
Version 1.0
Summer 2019


Table of Contents:

1. Introduction	2
2. Kit Contents	3
3. Storage and Stability	3
4. Additionally Required Materials	3
5. General Precautions	4
6. Specimen, storage and transport.....	4
7. Interfering substances	4
8. Internal control (IC).....	5
9. DNA isolation.....	5
10. PCR Protocol.....	6
11. Programming of the Rotor-Gene	6
12. Programming of StepOne.....	6
13. Programming other machines	7
14. Data Analysis: Rotor-Gene.....	7
15. Data Analysis: StepOne	9
16. Quantitation.....	10
17. Linear Range.....	12
18. Sensitivity.....	12

JC RQ kit is intended for the quantitative detection of JC virus DNA with Rotor-Gene or StepOne machines. This kit is for research use only.

1. Introduction

JC virus is a member of polyomaviruses genus with a double stranded DNA genome. BK virus belongs to this genus as well. Early infection occurs in early childhood and without clinical symptoms. Then JC virus establishes a lifelong infection usually in kidney, central nervous system or lymphocytes. Usually about 80% of adult population are infected with this virus without noticing it. However, JCV may lead to significant consequences in patients with immune deficiency or very weak immune system such as people with HIV/AIDs, lymphoma, multiple sclerosis (MS) and transplant patients. There, JCV reactivation may result on a serious brain infection called progressive multifocal leukoencephalopathy (PML). Therefore, it is crucial to monitor JCV titre in above mentioned target group.

JC RQ kit provides a ready-to-use Real-Time PCR system for detection and quantitation of JC virus DNA with Rotor-Gene or StepOne machines. Currently the Real-Time PCR provides the highest sensitivity and widest dynamic range among other methods. In this method, application of fluorescent probes allows detection of amplified product. Analysis of fluorescent kinetics also leads to quantification of the target sequence in the reaction without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product.

This kit also incorporates an Internal Control to identify possible PCR inhibition.

2. Kit Contents

The kit contains a manual, a CD with RotorGene and StepOne templates and following reagents:

Label	Content	Quantity
JC Mix	PCR Master mix*	360 μ l
JC S1	Standard 1: 100,000 copy/ μ l	150 μ l
JC S2	Standard 2: 10,000 copy/ μ l	150 μ l
JC S3	Standard 3: 1,000 copy/ μ l	150 μ l
JC S4	Standard 4: 100 copy/ μ l	150 μ l
JC IC	Internal Control	250 μ l
Water	PCR Grade Water	200 μ l

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

3. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws more than few times to prevent reduced sensitivity.

4. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Table top microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.7ml microtubes and PCR microtubes
- Disposable powder-free gloves

5. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the PCR mix is aliquoted into tubes and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working.**
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

6. Specimen, Storage and Transport

Peripheral blood, cerebrospinal fluid or urine should be collected in sterile condition in proper and sterile tubes. Samples can be shipped and stored at +4°C for couple of days or at -20°C for longer periods. Upon receipt plasma should be separated from whole blood and then stored as above mentioned.

7. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes should not be used. Samples of heparinized patients must not be used as well.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

8. Internal Control (IC)

In order to evaluate the possibility of DNA extraction failure and PCR inhibition and prevent false negative results, JC RQ kit contains internal control (IC). IC can be used during extraction process or simply added to JC Mix.

To monitor both DNA extraction and PCR reaction, internal control should be added to the mixture of lysis buffer and patient sample. Required volume of internal control is 10% of elution buffer. For instance if extracted DNA is eluted with 100ul, then 10ul of internal control should be added to the mixture of lysis buffer and patient sample. **Please note that internal control should not be added directly to the patient sample (i.e. before addition of lysis buffer) as it loses its efficiency.**

If IC is added to JC Mix, only PCR inhibition can be monitored. For this purpose, 1ul of IC reagent should be added to each reaction. For example, for 10 PCR reactions, 15ul of internal control should be added to 150ul of JC Mix before it is added to the tubes.

In a successful DNA extraction and PCR test, IC should generate a CT of 27-35 in Yellow/VIC Channel.

9. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend using either of:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)
- QIAamp UltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)

- QIAamp MiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

10. PCR Protocol

Thaw the reagents on ice completely followed by a brief mixing and a quick spin. Place required number of tubes on cold block. Consider one tube for each sample plus four for standards and one for the negative control.

Pipette 15ul of JC Mix directly to each tube followed by adding 10ul of standard or isolated DNA.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring.

11. Programming Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the CD provided with the kit and double click on JC template 0.2 or JC 0.1 depending on the tubes used. Program starts. Click on Start button (Green button on the top menu). On the pop up window, click start again and save program on desired location.

12. Programming StepOne

Open the StepOne software (V 2.*). On the Set Up menu click on Template and select the file on CD provided with the kit. Click on Plate Setup. One negative control, 4 standards and 10 samples are defined. You may change plate set up using right click options (copy, paste, clear). You may also add/remove samples or

change sample names on “Define Targets and Samples” menu. When finished, click on “Start Run” and save the experiment on desired location. Instrument will start shortly.

13. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 10 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. The JC Mix contains ROX with the final concentration of 300nM in the reaction.

14. Data Analysis: Rotor-Gene

Before analyzing results, make sure in the sample menu all the standards have been defined as “standard” and the relative concentrations have been entered. Patient samples should be defined as “unknown” and Negative control or no template control as “Negative Control” or “NTC” respectively.

Analyze data according to manufacturer recommendations. Perform quantitative analysis for **JC (Green channel)** and qualitative analysis for **Internal Control (Yellow channel)**. Briefly, click on analysis menu and then under Quantitation tab double click on cycling A. Green.

Close the pop up for Automatic and set the threshold at 0.1. Repeat the above for Cycling A. Yellow.

Figures 1 and 2 represent typical graphs for Rotor-Gene machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used.

In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.

Consider following points when analyzing:

- A sample is **Positive** if it is positive in Green channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results in the Cycling A. Green are valid.
- A sample is **Negative** if it is negative in Green channel while it is positive in Yellow channel with a sigmoid graph and CT of 28-35.
- Results are **inconclusive** and the test should be repeated if a sample is negative in both of Green and Yellow channels.


Fig 1. Typical JC graph in Green channel for Rotor-Gene


Fig2. Typical IC graph in Yellow channel for Rotor-Gene

15. Data Analysis: StepOne

Analyze data according to manufacturer recommendations. Briefly, click on Analyze and set the threshold for **JC/FAM** at 0.1 and at 0.05 for **IC/VIC**.

Figures 3 and 4 represent typical graphs for StepOne machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used.

In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.

Consider following points when analyzing:

- A sample is **Positive** if it is positive in FAM channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results are valid.
- A sample is **Negative** if it is negative in FAM channel while it is positive in VIC channel with a sigmoid graph and CT of 28-35.
- Results are **inconclusive** and the test should be repeated if a sample is negative in both of FAM and VIC channels.

16. Quantitation

The kit provides 4 quantitation standards with defined titers to generate a standard curve for quantification of samples viral load. Working with RotorGene machine, the standard curve from a previous run can also be imported for quantification of samples to the recent run. To do so, at least one standard must be used in the current run. Apparently using all four standards in each run will lead to more accurate results.

Quantitation standards are defined as copy/ μ l. To convert the result to copy/ml following equation should be used:

$$\text{Result (copy/ml)} = \frac{\text{Result (copy/}\mu\text{l)} \times \text{elution volume (}\mu\text{l)}}{\text{sample volume (ml)}}$$

“Sample volume” is the sample volume used for DNA isolation and “Elution volume” is the volume of buffer or water used to elute or dissolve isolated DNA.


Fig 3. Typical JC graph in FAM channel for StepOne


Fig 4. Typical IC graph in VIC channel for StepOne

17. Linear Range

The linear range of the kit was assessed with dilution series of the cloned target and showed to be linear in the range of 100,000,000 copy/ μ l to 50 copy/ μ l.

18. Sensitivity

The analytical detection limit of the kit was assessed with dilution series of the cloned target and showed a limit of detection equal to 5 copy/ μ l.