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HPV HR RQ3 Kit Manual

**For Real-Time PCR Detection of
14 High Risk and 2 Low Risk Human Papilloma
Virus (HPV) Types**

For use with Rotor-Gene or StepOne
Research use only

NG-WI-ASL-36-00
Version 1.1
Winter 2020


Table of Contents:

1. Introduction.....	2
2. Kit Contents.....	2
3. Storage and Stability.....	3
4. Additionally Required Materials.....	3
5. General Precautions.....	3
6. Specimen, storage and transport	4
7. DNA isolation.....	4
8. Internal control (IC).....	4
9. PCR Protocol.....	5
10. Programming Rotor-Gene	6
11. Programming StepOne	6
12. Programming Other Machines	6
13. Data Analysis: Rotor-Gene.....	7
14. Data Analysis: StepOne.....	11
15. Sensitivity.....	15

HPV HR RQ3 kit is intended for use with Rotor-Gene or StepOne machines and for the detecting DNA of 14 High Risk Types (16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68) and 2 Low Risk Types (6 and 11) of Human Papilloma virus. This kit is for research use only.

1. Introduction

Human Papilloma Virus (HPV) is a Non-enveloped double stranded DNA virus from papillomaviridae family. There are about 150 different types of HPV. HPV infection is considered the main cause of cervical cancer and has been linked to more than 99% of all cases. According to data linking HPV infection to cervical cancer, HPV types have been divided into three groups of High-Risk, Low-Risk and No-Risk. Currently 14 types of 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68 are considered High-Risk and 6 types considered Low-Risk including 6, 11, 42, 43, 44 and 53.

HPV HR RQ3 kit provides a ready-to-use Real-Time PCR system for detection of HPV DNA including 14 High Risk Types of 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68 and 2 Low Risk Types of 6 and 11 with RotorGene or StepOne machines. In this method, application of fluorescent dye labeled probes allows detection of amplified product without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product.

This kit also incorporates an *Internal Control* (IC) to evaluate DNA extraction quality or PCR inhibition.

2. Kit Contents

The kit contains a manual, a CD with Rotor-Gene and StepOne templates and following reagents:

Label	Content	Quantity
HPV HR1 Mix	PCR Master mix*	480 µl
HPV HR2 Mix	PCR Master mix*	480 µl
HPV HR3 Mix	PCR Master mix*	480 µl
HPV HR4 Mix	PCR Master mix*	480 µl
Pos Ctrl	Positive Control	100 µl
Neg Ctrl	Negative Control	100 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

3. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws more than few times to prevent reduced sensitivity.

4. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Table top microtube centrifuge
- Vortex Mixer
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.7ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

5. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place 0.2ml PCR tubes on crushed ice. Use cold blocks instead.

6. Specimen, Storage and Transport

Pap-smear or collection of cervicovaginal cells can be used for DNA extraction. Sample should be shipped at +4°C and upon receipt should be stored at +4°C for few hours or aliquoted and stored at -20°C for up to few weeks.

7. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure PCR Template Preparation Kit (Cat# 11796828001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)

8. Internal Control (IC)

To examine DNA extraction quality as well as the presence of PCR inhibitors and to prevent false negative results, primers and probe for an *Internal Control* (a housekeeping gene) is included in HPV HR1 Mix. With a proper sample, successful DNA extraction and PCR, Internal control should generate a CT of 20-30. Refer to analysis section for more details.

9. PCR Protocol

This kit contains 4 different PCR mixes, each for detection of 3 to 6 of HPV types. Therefore, each sample should be examined with all 4 mixes.

First, thaw the reagents on ice completely followed by a brief mixing and a quick spin. Place 4 series of tubes on cold block. Each series should include one tube for each sample plus three for positive control, negative control and no template control (NTC).

For example for examining 5 samples, 20 tubes for samples and 12 tubes for controls are required. There for a total of 32 tubes arranged in 4 series of 8 tubes each.

Pipette 20µl of HPV HR1 Mix to each of PCR tubes in the first series.

Pipette 20µl of HPV HR2 Mix to each of PCR tubes in the second series.

Pipette 20µl of HPV HR3 Mix to each of PCR tubes in the third series.

Pipette 20µl of HPV HR4 Mix to each of PCR tubes in the forth series.

Then add 5ul of extracted DNA, Pos Ctrl, Neg ctrl or water to each tube.

consider the first tube in each series for the first sample, second tube in each series for the second sample and continue respectively for the remaining samples or for the controls.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using RotorGene attach the locking ring.

10. Programming Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

CD provided in the kit. Open the HPV HR RQ3 folder then double click on HPV HR RQ3 0.2 or HPV HR RQ3 strip template depending on used tubes. Program starts. Click on Start button (Green button on the top menu). On the pop up window, click start again and save program on desired location.

11. Programming StepOne

Open the StepOne software (V 2.*). On the Set Up menu click on Template and select the HPV HR RQ3 file on CD provided with the kit. Click on Plate Setup. Controls and few samples are defined. You may change plate set up using right click options (copy, paste, clear). You may also add /remove samples or change sample name on "Define Targets and Samples" menu. When finished, click on "Start Run" and save the experiment on desired location. Instrument will start shortly.

12. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 10 min	1 cycle
2	95°C x 15 sec	45 cycles
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM, VIC and ROX dyes.

Please note that, HPV HR Mixes does not contain ROX dye as normalizer!

13. Data Analysis: Rotor-Gene

Analyze data according to manufacturer recommendations. Briefly, click on analysis menu and then under Quantitation tab double click on cycling A. Green. Close the pop up window for Automatic Threshold and set the threshold at 0.1.

Repeat the above for Cycling A. Yellow and Orange. To interpret the results, please, note the following table. It shows each of the 4 mixes, detects which of the HPV in each channel.

	Green	Yellow	Orange
HR1 Mix	16	6/11	Albumin
HR2 Mix	18	39	59
HR3 Mix	31/33	52/58	56/66
HR4 Mix	35/45	51	68

According to the above table, results can be interpreted as below:

Note that a sample is considered Positive only if it has a log phase and sigmoid graph with a CT less than 40. In the absence of log phase and a sigmoid graph and with a CT of above 40, sample can not be considered as Positive.

Please note that, for a valid result, each sample should be positive in Orange channel with CT of 20-30. When a sample is Negative in Orange channel or has a CT above 30, results are not valid and test should be repeated.

Negative Result: A sample is **Negative** if it is **only** positive with **HR1 Mix** in **Orange** channel and remains negative for all other mixes-channels combinations.

Positive Result: If a sample is positive with HR1 Mix in **Orange** channel with CT of 20-30, then:

- It is **Positive for HPV16** if it is positive in **Green** channel with **HR1 Mix**.
- It is **Positive for HPV6 or 11** if it is positive with CT<35 in **Yellow** channel with **HR1 Mix**.
- It is **Positive for HPV18** if it is positive in **Green** channel with **HR2 Mix**.
- It is **Positive for HPV39** if it is positive in **Yellow** channel with **HR2 Mix**.
- It is **Positive for HPV59** if it is positive in **Orange** channel with **HR2 Mix**.
- It is **Positive for HPV31 or 33** if it is positive in **Green** channel with **HR3 Mix**.
- It is **Positive for HPV52 or 58** if it is positive in **Yellow** channel with **HR3 Mix**.
- It is **Positive for HPV56/66** if it is positive in **Orange** channel with **HR3 Mix**.

- It is **Positive for HPV35 or 45** if it is positive in **Green** channel with **HR4 Mix**.
- It is **Positive for HPV51** if it is positive in **Yellow** channel with **HR4 Mix**.
- It is **Positive for HPV68** if it is positive in **Orange** channel with **HR4 Mix**.

Above results are summarized in following table:

	Green	Yellow	Orange	Result
HR1 Mix	+	-	IC	Pos: HPV16
HR1 Mix	-	+ CT<35	IC	Pos: HPV6 or 11
HR2 Mix	+	-	-	Pos: HPV 18
HR2 Mix	-	+	-	Pos HPV39
HR2 Mix	-	-	+	Pos: HPV59
HR3 Mix	+	-	-	Pos: HPV31 or 33
HR3 Mix	-	+	-	Pos: HPV52 or 58
HR3 Mix	-	-	+	Pos: HPV56 or 66
HR4 Mix	+	-	-	Pos: HPV35 or 45
HR4 Mix	-	+	-	Pos: HPV51
HR4 Mix	-	-	+	Pos: HPV68

Figures 1, 2 and 3 represent typical graphs for Rotor-Gene machine.


Fig 1. Typical HPV HR graph in Green channel for Rotor-Gene


Fig 2. Typical HPV HR graph in Yellow channel for Rotor-Gene


Fig 3. Typical HPV HR graph in Orange channel for Rotor-Gene

14. Data Analysis: StepOne

Analyze data according to manufacturer recommendations. Briefly, click on Analyze and set the threshold for FAM and VIC at 1000 and at 500 for ROX. To interpret the results, please, note the following table. It shows each of the 4 mixes, detects which of the HPV in each channel.

	FAM	VIC	ROX
HR1 Mix	16	6/11	Albumin
HR2 Mix	18	39	59
HR3 Mix	31/33	52/58	56/66
HR4 Mix	35/45	51	68

According to the above table, results can be interpreted as below:

Note that a sample is considered Positive only if it has a log phase and sigmoid graph with a CT less than 40. In the absence of log phase and a sigmoid graph and with a CT of above 40, sample can not be considered as Positive.

Please note that, for a valid result, each sample should be positive in ROX channel with CT of 20-30. When a sample is Negative in ROX channel or has a CT above 30, results are not valid and test should be repeated.

Negative Result: A sample is **Negative** if it is **only** positive with **HR1 Mix** in **ROX** channel and remains negative for all other mixes-channels combinations.

Positive Result: If a sample is positive with HR1 Mix in **Orange** channel with CT of 20-30, then:

- It is **Positive for HPV16** if it is positive in **FAM** channel with **HR1 Mix**.
- It is **Positive for HPV6 or 11** if it is positive with CT<35 in **VIC** channel with **HR1 Mix**.
- It is **Positive for HPV18** if it is positive in **FAM** channel with **HR2 Mix**.
- It is **Positive for HPV39** if it is positive in **VIC** channel with **HR2 Mix**.
- It is **Positive for HPV59** if it is positive in **ROX** channel with **HR2 Mix**.
- It is **Positive for HPV31 or 33** if it is positive in **FAM** channel with **HR3 Mix**.
- It is **Positive for HPV52 or 58** if it is positive in **VIC** channel with **HR3 Mix**.
- It is **Positive for HPV56/66** if it is positive in **ROX** channel with **HR3 Mix**.

- It is **Positive for HPV35 or 45** if it is positive in **FAM** channel with **HR4**.
- It is **Positive for HPV51** if it is positive in **VIC** channel with **HR4 Mix**.
- It is **Positive for HPV68** if it is positive in **ROX** channel with **HR4 Mix**.

Above results are summarized in following table:

	FAM	VIC	ROX	Result
HR1 Mix	+	-	IC	Pos: HPV16
HR1 Mix	-	+ CT<35	IC	Pos: HPV6 or 11
HR2 Mix	+	-	-	Pos: HPV 18
HR2 Mix	-	+	-	Pos HPV39
HR2 Mix	-	-	+	Pos: HPV59
HR3 Mix	+	-	-	Pos: HPV31 or 33
HR3 Mix	-	+	-	Pos: HPV52 or 58
HR3 Mix	-	-	+	Pos: HPV56 or 66
HR4 Mix	+	-	-	Pos: HPV35 or 45
HR4 Mix	-	+	-	Pos: HPV51
HR4 Mix	-	-	+	Pos: HPV68

Figures 4, 5 and 6 represent typical graphs for StepOne machine.


Fig 4. Typical HPV HR graph in FAM channel for StepOne


Fig 5. Typical HPV HR graph in VIC channel for StepOne


Fig 6. Typical HPV HR graph in ROX channel for StepOne

15. Sensitivity

The analytical detection limit of the kit was assessed with dilution series of the cloned HPV 16 target and the detection limit was determined as 50 copy/ μ l.