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
# **HCV Genotype TM**

## **Kit Manual**

**For Genotyping of Hepatitis C Virus  
(HCV) with Real-Time RT-PCR**

For use with StepOne  
Research use only

NG-WI-ASL-04-200  
Version 2.0  
March 2017

## Table of Contents:

1. Introduction.....	2
2. Kit Contents.....	3
3. Storage and Stability.....	3
4. General Precautions.....	3
5. Additionally Required Materials.....	4
6. Specimen, Storage and Transport.....	4
7. Interfering Substances.....	5
8. RNA Isolation.....	5
9. RT-PCR Protocol.....	5
10. Programming of StepOne.....	6
11. Programming Other Machines.....	6
12. Data Analysis: StepOne.....	7
13. Sensitivity.....	10

**HCV Gen TM** kit is intended for the genotyping of HCV types 1, 2 and 3 using RNA extracted from plasma. It is designed for use with StepOne instrument. This kit is for research use only.

## 1. Introduction

Hepatitis C virus (HCV) is the major cause of cirrhosis and hepatocellular carcinoma. So far, 6 major genotypes and more than 80 subtypes have been identified.

Based on current clinical trials, pre-treatment viral load and the HCV genotype are two best indicators for medication choice as well as patient response to therapy. Therefore, once patient tested positive for HCV, both viral load and genotype should be determined. Moreover since HCV genotypes have significant differences in their global distribution and prevalence, genotyping can be a useful tool for epidemiologic studies.

The gold standard method for HCV genotyping is sequencing and phylogenetic studies. Since this method is not feasible for routine practice in a clinical laboratory, other methods have been used for this purpose which most of them are PCR-based.

HCV Genotype TM kit provides a ready-to-use One-Step Real-Time RT-PCR system for genotyping of HCV types 1, 2 and 3. This method application of fluorescent probes provides data about reaction kinetics. Analysis of fluorescent kinetics leads to detection and genotyping of the target sequence in the reaction without requiring post-amplification analysis therefore, reducing the possibility of contamination with the PCR product. This kit is designed to be used with StepOne instrument. This kit does not detect other HCV types including types 4, 5 and 6.

## 2. Kit Contents

The kit contains a manual, a CD with StepOne template and following reagents:

Label	Content	Quantity
HCV Gen Mix TM	RT-PCR Mix*	300µl
HCV-1	HCV Type 1 control	150µl
HCV-2	HCV Type 2 control	150µl
HCV-3	HCV Type 3 control	150µl
Water	PCR Grade Water	200µl

\* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

## 3. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws especially for HCV Gen Mix more than few times to prevent reduced sensitivity.

## 4. General Precautions

In order to prevent false results, always pay attention to the following points:

- Treat all samples as potentially infectious.
- Within the pre-PCR work area assign three separate spaces for: a) Sample storage and extraction;; b) Reagent preparation where the master-mix is aliquoted into tubes; and c) Reaction preparation area for addition of extracted RNA to the tubes.

- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on crushed ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place 0.2ml PCR tubes on crushed ice. Use cold blocks instead.

## 5. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Table top microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- RNA extraction kit
- Nuclease free 1.7ml microtubes and PCR microtubes
- Disposable powder-free gloves

*Note! This kit does not require cDNA synthesis reagents. They are already included in the Mix.*

## 6. Specimen, Storage and Transport

We recommend EDTA or citrate plasma for HCV detection and typing. Peripheral blood should be collected in sterile condition in proper and sterile tubes. Whole blood or plasma should be shipped at +4°C. Upon receipt, plasma should be separated from whole blood and can be stored at +4°C for few days or aliquoted and stored at -20°C for up to few weeks.

## 7. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes should not be used. Samples of heparinized patients must not be used as well. Elevated levels of bilirubin ( $\leq 4.5$  mg/dl) and lipids ( $\leq 1000$  mg/dl) and hemolytic samples do not influence the extraction and PCR.

## 8. RNA Isolation

RNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp Viral RNA Mini Kit (Cat. no. 52904, Qiagen GmbH, Hilden, Germany)
- QIAamp UltraSens® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAamp MiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

## 9. RT-PCR Protocol

Thaw the reagents on crushed ice completely followed by a brief mixing and a quick spin. Place required number of tubes on cold block. Consider one tube for each sample plus four for positives and negative control.

**Pipette 12.5µl of HCV Gen TM Mix directly to each tube followed by adding 12.5µl of isolated RNA or control.**

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

*Note: Working with StepOne instrument, spin tubes briefly before loading on the block.*

### 10. Programming StepOne

Open the StepOne software (V 2.\*). On the Set Up menu click on Template and select the file on CD provided with the kit. Click on Plate Setup. Three positive controls for HCV types, one negative control and a few samples are defined. You may change plate set up using right click options (copy, past, clear). You may also add/remove samples or change sample name on “Define Targets and Samples” menu. When finished, click on “Start Run” and save the experiment on desired location. Instrument will start shortly.

*Note! ROX should not be selected as reference dye in this test.*

### 11. Programming Other Machines

**HCV Gen TM** kit is intended for use with StepOne machine. If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	<b>50°C x 30 min</b>	1
2	<b>95°C x 15 min</b>	1
3	<b>95°C x 15 sec</b>	10
	<b>63°C x 15 sec</b>	
	<b>72°C x 15 sec</b>	
4	<b>95°C x 15 sec</b>	50
	<b>62°C x 15 sec</b>	
	<b>72°C x 15 sec</b>	


Fluorescence should be collected at 62°C for FAM, VIC and ROX during the final 50 cycles (section 4 of the profile).

*Note! ROX should not be selected as reference dye in this test.*


## 12. Data Analysis: StepOne

Analyze data according to manufacturer recommendations. Briefly set threshold to 1000 for all three channels of FAM, VIC and ROX.


For quick and easy interpretation of results, in analysis menu select the sample to view its graph. In the graph submenu named as “Target”, select the item “all”. Now, the lowest CT shows the sample type. Note that a signal in **FAM channel is due to Type 1**, in **VIC channel due to Type 3** and in **ROX channel due to Type 2**. Figures 1, 2 and 3 represent typical graphs of 3 genotypes for StepOne machine.

**Figure 1:** HCV Type 1 graph in FAM, VIC & ROX Channels

**Figure 2:** HCV Type 2 graph in FAM, VIC & ROX Channels

**Figure 3:** HCV Type 3 graph in FAM, VIC & ROX Channels

**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used.**

**In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.**

For more details consider the following points:

- Sample contains **HCV Type 1**, if it is positive in **FAM** channel. It may react in **VIC** channel with a CT of about 2 units or more above **FAM** CT.
- Sample contains **HCV Type 3**, if it is positive in **VIC** channel. It may react in **ROX** channel with a CT of about 6 units or more above **VIC** CT.
- Sample contains **HCV Type 2**, if it is positive in **ROX** channel.
- Sample contains both **HCV Type 1** and **HCV Type 3**, if it is positive in both of **FAM** and **VIC** channels while the CT of **VIC** channel is less than 2 units apart from **FAM** CT.

If sample is negative in all three channels of FAM, VIC and ROX, there are following possibilities:

- Sample is negative for HCV.
- Sample is positive for HCV genotypes 1, 2 and 3, but the viral load is lower than the test sensitivity.
- Sample is positive for other HCV genotypes including 4, 5 or 6.

*Note: according to published studies, type 1 is the most prevalent type in Iran with type 3 as the second most prevalent type. Type 2 has also been reported. Types 4, 5 and 6 have not been reported from Iran yet.*

### **13. Sensitivity**

The analytical detection limit of the kit was assessed with dilution series of the cloned target (partial type 1 genome) and showed a limit of detection equal to 6 IU/ $\mu$ l.