

BIO BASIC INC.
20 KONRAD CRES, MARKHAM ONTARIO L3R 8T4 CANADA
TEL: (905) 474 4493, (800) 313 7224 FAX: (905) 474 5794

MATERIAL SAFETY DATA SHEET
REQUIRED UNDER SAFETY AND HEALTH REGULATION FOR SHIP REPAIRING

MANUFACTURER'S NAME: BIO BASIC INC.
TEL: (905) 474 4493
FAX: (905) 474 5794
FOR CHEMICAL EMERGENCY: (416) 995 9730
ADDRESS: 20 KONRAD CRES, MARKHAM ONTARIO L3R 8T4 CANADA

DATE UPDATED: OCTOBER 28, 2014

SECTION 1. ----- CHEMICAL IDENTIFICATION -----

Product Name Amoxicillin trihydrate
Product Code(s) A762203
Recommended Use For Further Manufacturing Use Only
Not for Human or Animal Drug Use

SECTION 2. ----- HAZARDS IDENTIFICATION -----

Emergency Overview

WHMIS Classification

Not WHMIS controlled.

GHS Classification

Respiratory sensitisation (Category 1)
Skin sensitisation (Category 1)

GHS Label elements, including precautionary statements

Pictogram



Signal word Danger

Hazard statement(s)

H317 May cause an allergic skin reaction.

H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Precautionary statement(s)

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

P280 Wear protective gloves.

P342 + P311 If experiencing respiratory symptoms: Call a POISON CENTER or doctor/ physician.

HMIS Classification

Health hazard: 0

Flammability: 0

Physical hazards: 0

Potential Health Effects

Inhalation May be harmful if inhaled. May cause respiratory tract irritation.

Skin May be harmful if absorbed through skin. May cause skin irritation.

Eyes May cause eye irritation.

Ingestion May be harmful if swallowed.

SECTION 3. - - - - - COMPOSITION/INFORMATION ON INGREDIENTS - - - - -

Chemical Name	EC No.	REACH Reg. No.	CAS-No	Weight %	Classification
Amoxicillin trihydrate	-	Not available	61336-70-7	95-100	-

SECTION 4. - - - - - FIRST-AID MEASURES - - - - -**General advice**

Consult a physician. Show this safety data sheet to the doctor in attendance.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

In case of skin contact

Wash off with soap and plenty of water. Consult a physician.

In case of eye contact

Flush eyes with water as a precaution.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

SECTION 5. - - - - - FIRE FIGHTING MEASURES - - - - -**Conditions of flammability**

Not flammable or combustible.

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special protective equipment for firefighters

Wear self-contained breathing apparatus for firefighting if necessary.

Hazardous combustion products

Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NO_x), Sulphur oxides

Explosion data - sensitivity to mechanical impact

No data available

Explosion data - sensitivity to static discharge

No data available

SECTION 6. - - - - - ACCIDENTAL RELEASE MEASURES - - - - -**Personal precautions**

Use personal protective equipment. Avoid dust formation. Avoid breathing vapours, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas. Avoid breathing dust.

Environmental precautions

Do not let product enter drains.

Methods and materials for containment and cleaning up

Pick up and arrange disposal without creating dust. Sweep up and shovel. Keep in suitable, closed containers for disposal.

SECTION 7. - - - - - HANDLING AND STORAGE - - - - -**Precautions for safe handling**

Avoid contact with skin and eyes. Avoid formation of dust and aerosols.

Provide appropriate exhaust ventilation at places where dust is formed. Normal measures for preventive fire protection.

Conditions for safe storage

Keep container tightly closed in a dry and well-ventilated place.

Recommended storage temperature 2 - 8 °C

SECTION 8. - - - - - EXPOSURE CONTROLS/PERSONAL PROTECTION- - - - -**Personal protective equipment****Respiratory protection**

For nuisance exposures use type P95 (US) or type P1 (EU EN 143) particle respirator. For higher level protection use type OV/AG/P99 (US) or type ABEK-P2 (EU EN 143) respirator cartridges. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Eye protection

Face shield and safety glasses Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Skin and body protection

Complete suit protecting against chemicals, The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Hygiene measures

Handle in accordance with good industrial hygiene and safety practice. Wash hands before breaks and at the end of workday.

Specific engineering controls

Use mechanical exhaust or laboratory fumehood to avoid exposure.

SECTION 9. - - - - - PHYSICAL AND CHEMICAL PROPERTIES - - - - -**Appearance**

Form	solid
Colour	No data available

Safety data

pH	No data available
Melting point/freezing point	No data available
Boiling point	No data available
Flash point	No data available
Ignition temperature	No data available
Auto-ignition temperature	No data available
Lower explosion limit	No data available
Upper explosion limit	No data available
Vapour pressure	No data available
Density	No data available
Water solubility	No data available
Partition coefficient: n-octanol/water	No data available

Relative vapour density	No data available
Odour	No data available
Odour Threshold	No data available
Evaporation rate	No data available

SECTION 10. -----STABILITY AND REACTIVITY -----

Chemical stability

Stable under recommended storage conditions.

Possibility of hazardous reactions

No data available

Conditions to avoid

No data available

Materials to avoid

Strong oxidizing agents

Hazardous decomposition products

Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NOx), Sulphur oxides
Other decomposition products - No data available

SECTION 11. ----- TOXICOLOGICAL INFORMATION -----

Acute toxicity

Oral LD50

LD50 Oral - Rat - > 15,000 mg/kg

Inhalation LC50

No data available

Dermal LD50

No data available

Other information on acute toxicity

No data available

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

May cause allergic respiratory reaction. May cause allergic skin reaction. Shows cross-allergenicity with penicillin.

May cause allergic respiratory and skin reactions

Germ cell mutagenicity

No data available

Carcinogenicity

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

Reproductive toxicity

No data available

Teratogenicity

Developmental Toxicity - Rat - Oral

Specific Developmental Abnormalities: Musculoskeletal system.

Developmental Toxicity - Mouse - Oral

Effects on Embryo or Fetus: Fetal death. Specific Developmental Abnormalities: Musculoskeletal system.

No data available

Specific target organ toxicity - single exposure (Globally Harmonized System)

No data available

Specific target organ toxicity - repeated exposure (Globally Harmonized System)

No data available

Aspiration hazard

No data available

Potential health effects

Inhalation	May be harmful if inhaled. May cause respiratory tract irritation.
Ingestion	May be harmful if swallowed.
Skin	May be harmful if absorbed through skin. May cause skin irritation.
Eyes	May cause eye irritation.

Synergistic effects

No data available

Additional Information

RTECS: XH8310000

SECTION 12. ----- ECOLOGICAL INFORMATION -----

Toxicity

No data available

Persistence and degradability

No data available

Bioaccumulative potential

No data available

Mobility in soil

No data available

PBT and vPvB assessment

No data available

Other adverse effects

No data available

SECTION 13. ----- DISPOSAL CONSIDERATIONS -----

Product

Offer surplus and non-recyclable solutions to a licensed disposal company.

Contaminated packaging

Dispose of as unused product.

14. TRANSPORT INFORMATION

DOT (US)

Not dangerous goods

IMDG

Not dangerous goods

IATA

Not dangerous goods

SECTION 15. ----- REGULATORY INFORMATION -----**WHMIS Classification**

Not WHMIS controlled.

This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations and the MSDS contains all the information required by the Controlled Products Regulations.

SECTION 16. ----- OTHER INFORMATION -----

Issuing Date	13-Aug-2009
Revision Date	28-Oct-2014
Revision Note	No information available.
Recommended Restrictions	No information available

Disclaimer

The information provided on this MSDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

End of MSDS