

BIO BASIC INC.
20 KONRAD CRES, MARKHAM ONTARIO L3R 8T4 CANADA
TEL: (905) 474 4493, (800) 313 7224 FAX: (905) 474 5794

MATERIAL SAFETY DATA SHEET
REQUIRED UNDER SAFETY AND HEALTH REGULATION FOR SHIP REPAIRING

MANUFACTURER'S NAME: BIO BASIC INC.
TEL: (905) 474 4493
FAX: (905) 474 5794
FOR CHEMICAL EMERGENCY: (416) 995 9730
ADDRESS: 20 KONRAD CRES, MARKHAM ONTARIO L3R 8T4 CANADA

DATE UPDATED: JANUARY 08, 2015

SECTION 1. ----- CHEMICAL IDENTIFICATION-----

Product Name Ammonium iron(II) sulfate hexahydrate
Product Code(s) AB0065
Recommended Use For Laboratory Use Only
Not for Human or Animal Drug Use

SECTION 2. ----- HAZARDS IDENTIFICATION -----

Emergency Overview

WHMIS Classification

D2B Toxic Material Causing Other Toxic Effects Moderate skin irritant
Moderate eye irritant

GHS Classification

Acute toxicity, Oral (Category 5)
Skin irritation (Category 2)
Eye irritation (Category 2A)
Specific target organ toxicity - single exposure (Category 3)

GHS Label elements, including precautionary statements

Pictogram



Signal word Warning

Hazard statement(s)

H303 May be harmful if swallowed.
H315 Causes skin irritation.
H319 Causes serious eye irritation.
H335 May cause respiratory irritation.

Precautionary statement(s)

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

HMIS Classification

Health hazard: 2
Flammability: 0
Physical hazards: 0

Potential Health Effects

Inhalation May be harmful if inhaled. Causes respiratory tract irritation.
Skin May be harmful if absorbed through skin. Causes skin irritation.
Eyes Causes eye irritation.

Ingestion

May be harmful if swallowed.

SECTION 3. - - - - - COMPOSITION/INFORMATION ON INGREDIENTS - - - - -

Synonyms : Ammonium iron(II) sulfate
Mohr's salt
Ammonium ferrous sulfate hexahydrate

Formula : $HgFeN_2O_8S_2 \cdot 6H_2O$
Molecular Weight : 392.14 g/mol

CAS-No.	EC-No.	Index-No.	Concentration
Ammonium ferrous sulfate hexahydrate			
7783-85-9	233-151-8	-	<=100%

SECTION 4. - - - - - FIRST-AID MEASURES - - - - -

General advice

Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

In case of skin contact

Wash off with soap and plenty of water. Consult a physician.

In case of eye contact

Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

SECTION 5. - - - - - FIRE FIGHTING MEASURES - - - - -

Conditions of flammability

Not flammable or combustible.

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special protective equipment for firefighters

Wear self contained breathing apparatus for fire fighting if necessary.

Hazardous combustion products

Hazardous decomposition products formed under fire conditions. - nitrogen oxides (NOx), Sulphur oxides, Iron oxides

Explosion data - sensitivity to mechanical

impact no data available

Explosion data - sensitivity to static

discharge no data available

SECTION 6. - - - - - ACCIDENTAL RELEASE MEASURES - - - - -

Personal precautions

Use personal protective equipment. Avoid dust formation. Avoid breathing vapours, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas. Avoid breathing dust.

Environmental precautions

Do not let product enter drains.

Methods and materials for containment and cleaning up

Pick up and arrange disposal without creating dust. Sweep up and shovel. Keep in suitable, closed containers for disposal.

SECTION 7. ----- HANDLING AND STORAGE-----

Precautions for safe handling

Avoid contact with skin and eyes. Avoid formation of dust and aerosols. Provide appropriate exhaust ventilation at places where dust is formed.

Conditions for safe storage

Keep container tightly closed in a dry and well-ventilated place.

Air and light sensitive. Store under inert gas.

SECTION 8. ----- EXPOSURE CONTROLS/PERSONAL PROTECTION-----

Components with workplace control parameters

Components	CAS-No.	Value	Control parameters	Basis
Ammonium ferrous sulfate hexahydrate	7783-85-9	TWA	1 mg/m3	Canada. Alberta, Occupational Health and Safety Code (table 2: OEL)
Remarks	Occupational exposure limit is based on irritation effects and its adjustment to compensate for unusual work schedules is not required			
		TWAEV	1 mg/m3	Québec. Regulation respecting occupational health and safety, Schedule 1, Part 1: Permissible exposure values for airborne contaminants
		TWA	1 mg/m3	Canada. British Columbia OEL
		STEL	2 mg/m3	Canada. British Columbia OEL
		TWA	1 mg/m3	USA. ACGIH Threshold Limit Values (TLV)

Personal protective equipment

Respiratory protection

For nuisance exposures use type P95 (US) or type P1 (EU EN 143) particle respirator. For higher level protection use type OV/AG/P99 (US) or type ABEK-P2 (EU EN 143) respirator cartridges. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Eye protection

Safety glasses with side-shields conforming to EN166 Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Skin and body protection

impervious clothing, The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Hygiene measures

Handle in accordance with good industrial hygiene and safety practice. Wash hands before breaks and at the end of workday.

Specific engineering controls

Use mechanical exhaust or laboratory fumehood to avoid exposure.

SECTION 9. - - - - - PHYSICAL AND CHEMICAL PROPERTIES - - - - -

Appearance

Form	solid
Colour	no data available

Safety data

pH	no data available
Melting point/freezing point	Melting point/range: 100 °C (212 °F) - dec.
Boiling point	no data available
Flash point	not applicable
Ignition temperature	no data available
Auto-ignition temperature	no data available
Lower explosion limit	no data available
Upper explosion limit	no data available
Vapour pressure	no data available
Density	no data available
Water solubility	no data available
Partition coefficient: n-octanol/water	no data available
Relative vapour density	no data available
Odour	no data available
Odour Threshold	no data available
Evaporation rate	no data available

SECTION 10. - - - - - STABILITY AND REACTIVITY - - - - -

Chemical stability

Stable under recommended storage conditions.

Possibility of hazardous reactions no data available

Conditions to avoid

Light. Air

Materials to avoid

Strong acids, Strong oxidizing agents

Hazardous decomposition products

Hazardous decomposition products formed under fire conditions. - nitrogen oxides (NOx), Sulphur oxides, Iron oxides Other decomposition products - no data available

SECTION 11. - - - - - TOXICOLOGICAL INFORMATION - - - - -

Acute toxicity

Oral LD50

no data available

Inhalation LC50

no data available

Dermal LD50

no data
available

Other information on acute toxicity no data available

Skin corrosion/irritation

no data available

Serious eye damage/eye irritation

no data available

Respiratory or skin sensitisation

no data available

Germ cell mutagenicity

no data available

Carcinogenicity

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

Reproductive toxicity

no data available

Teratogenicity

no data available

Specific target organ toxicity - single exposure (Globally Harmonized System)

Inhalation - May cause respiratory irritation.

Specific target organ toxicity - repeated exposure (Globally Harmonized System) no data available

Aspiration hazard

no data
available

Potential health effects

- Inhalation** May be harmful if inhaled. Causes respiratory tract irritation.
- Ingestion** May be harmful if swallowed.
- Skin** May be harmful if absorbed through skin. Causes skin irritation.
- Eyes** Causes eye irritation.

Signs and Symptoms of Exposure

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

Synergistic effects

no data
available

Additional Information

RTECS: BR6500000

SECTION 12. ----- ECOLOGICAL INFORMATION -----

Toxicity

no data available

Persistence and degradability

no data
available

Bioaccumulative potential no data available

Mobility in soil
no data available

PBT and vPvB assessment no data available

Other adverse effects
no data available

SECTION 13. ----- DISPOSAL CONSIDERATIONS -----

Product

Offer surplus and non-recyclable solutions to a licensed disposal company. Contact a licensed professional waste disposal service to dispose of this material.

Contaminated packaging

Dispose of as unused product.

SECTION 14. ----- TRANSPORT INFORMATION -----

DOT (US)

UN number: 3077 Class: 9 Packing group: III
Proper shipping name: Environmentally hazardous substance, solid, n.o.s. (Ammonium ferrous sulfate hexahydrate) Reportable Quantity (RQ): 1000 lbs
Marine pollutant: No
Poison Inhalation Hazard: No

IMDG

Not dangerous goods

IATA

Not dangerous goods

SECTION 15. ----- REGULATORY INFORMATION -----

WHMIS Classification

D2B	Toxic Material Causing Other Toxic Effects	Moderate skin irritant Moderate eye irritant
-----	--	---

This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations and the MSDS contains all the information required by the Controlled Products Regulations.

SECTION 16. ----- OTHER INFORMATION -----

Issuing Date	13-Aug-2009
Revision Date	08-Jan-2015
Revision Note	No information available.
Recommended Restrictions	No information available

Disclaimer

The information provided on this MSDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

End of MSDS