

BIO BASIC INC.

160 Torbay Road, Markham Ontario. L3R 1G6 Canada
Tel: (905) 474 4493, (800) 313 7224 Fax: (905) 474 5794
MATERIAL SAFETY DATA SHEET

REQUIRED UNDER SAFETY AND HEALTH REGULATION FOR SHIP REPAIRING

MANUFACTURER'S NAME: BIO BASIC INC.
TEL: (905) 474 5794
FAX: (416) 995 9730
FOR CHEMICAL EMERGENCY: (416) 995 9730
ADDRESS: 160 Torbay Road, Markham Ontario. L3R 1G6 Canada

Dec. 1st, 2009

SECTION 1. ----- CHEMICAL IDENTIFICATION -----

CODE: BB0072

NAME: BCIP, P-TOLUIDINE SALT

SECTION 2. ----- COMPOSITION/INFORMATION ON INGREDIENTS -----

CHEMICAL NAME: CAS#: RTECS #:

BCIP 6578-06-9 NOT RATED

FOR MORE INFORMATION SEE COMPLETE RTECS ENTRY.

EC NO: 229-506-1

SECTION 3. ----- HAZARDS IDENTIFICATION -----

LABEL PRECAUTIONARY STATEMENTS

IRRITANT

IRRITATING TO EYES, RESPIRATORY SYSTEM AND SKIN.

IN CASE OF CONTACT WITH EYES, RINSE IMMEDIATELY WITH PLENTY OF WATER AND SEEK MEDICAL ADVICE.

WEAR SUITABLE PROTECTIVE CLOTHING.

SECTION 4. ----- FIRST-AID MEASURES -----

IF SWALLOWED, WASH OUT MOUTH WITH WATER PROVIDED PERSON IS CONSCIOUS.
CALL A PHYSICIAN.

IN CASE OF SKIN CONTACT, FLUSH WITH COPIOUS AMOUNTS OF WATER FOR AT LEAST 15 MINUTES. REMOVE CONTAMINATED CLOTHING AND SHOES. CALL A PHYSICIAN.

IF INHALED, REMOVE TO FRESH AIR. IF BREATHING BECOMES DIFFICULT, CALL A PHYSICIAN.

IN CASE OF CONTACT WITH EYES, FLUSH WITH COPIOUS AMOUNTS OF WATER FOR AT LEAST 15 MINUTES. ASSURE ADEQUATE FLUSHING BY SEPARATING THE EYELIDS WITH FINGERS. CALL A PHYSICIAN.

SECTION 5. ----- FIRE FIGHTING MEASURES -----

EXTINGUISHING MEDIA

CARBON DIOXIDE, DRY CHEMICAL POWDER OR APPROPRIATE FOAM.

SPECIAL FIREFIGHTING PROCEDURES

WEAR SELF-CONTAINED BREATHING APPARATUS AND PROTECTIVE CLOTHING TO PREVENT CONTACT WITH SKIN AND EYES.

SECTION 6. ----- ACCIDENTAL RELEASE MEASURES -----

WEAR PROTECTIVE EQUIPMENT.

SWEEP UP, PLACE IN A BAG AND HOLD FOR WASTE DISPOSAL.

AVOID RAISING DUST.

VENTILATE AREA AND WASH SPILL SITE AFTER MATERIAL PICKUP IS COMPLETE.

SECTION 7. ----- HANDLING AND STORAGE-----

REFER TO SECTION 8.

SECTION 8. ----- EXPOSURE CONTROLS/PERSONAL PROTECTION-----

WEAR APPROPRIATE NIOSH/MSHA-APPROVED RESPIRATOR, CHEMICAL-RESISTANT GLOVES, SAFETY GOGGLES, OTHER PROTECTIVE CLOTHING.

MECHANICAL EXHAUST REQUIRED.

HYGROSCOPIC

SECTION 9. ----- PHYSICAL AND CHEMICAL PROPERTIES -----

APPEARANCE

POWDER

SECTION 10. ----- STABILITY AND REACTIVITY -----

STABILITY

STABLE.

HAZARDOUS COMBUSTION OR DECOMPOSITION PRODUCTS

TOXIC FUMES OF:

CARBON MONOXIDE, CARBON DIOXIDE

NITROGEN OXIDES

HYDROGEN CHLORIDE GAS

HYDROGEN BROMIDE GAS

HAZARDOUS POLYMERIZATION

WILL NOT OCCUR.

INCOMPATIBILITIES:

MOISTURE

SECTION 11. ----- TOXICOLOGICAL INFORMATION -----

ACUTE EFFECTS

CAUSES EYE AND SKIN IRRITATION.

MATERIAL IS IRRITATING TO MUCOUS MEMBRANES AND UPPER RESPIRATORY TRACT.

MAY BE HARMFUL IF ABSORBED THROUGH THE SKIN, INHALED OR SWALLOWED.

THE TOXICOLOGICAL PROPERTIES HAVE NOT BEEN THOROUGHLY

INVESTIGATED.

SECTION 12. ----- ECOLOGICAL INFORMATION -----

DATA NOT YET AVAILABLE.

SECTION 13. ----- DISPOSAL CONSIDERATIONS -----

OBSERVE ALL FEDERAL, STATE AND LOCAL ENVIRONMENTAL REGULATIONS.

SECTION 14. ----- TRANSPORT INFORMATION -----

NOT REGULATED.

SECTION 15. ----- REGULATORY INFORMATION -----

EUROPEAN INFORMATION

IRRITANT

R: 36/37/38

IRRITATING TO EYES, RESPIRATORY SYSTEM AND SKIN.

S: 36/37/39

WEAR SUITABLE PROTECTIVE CLOTHING, GLOVES AND EYE/FACE PROTECTION

SECTION 16. ----- OTHER INFORMATION-----

THE ABOVE INFORMATION IS BELIEVED TO BE CORRECT BUT DOES NOT PURPORT TO BE ALL INCLUSIVE AND SHALL BE USED ONLY AS A GUIDE. AMRESKO INC. SHALL NOT BE HELD LIABLE FOR ANY DAMAGE RESULTING FROM HANDLING OR FROM CONTACT WITH THE ABOVE PRODUCT. SEE REVERSE SIDE OF INVOICE OR PACKING SLIP FOR ADDITIONAL TERMS AND CONDITIONS OF SALE.