

SRS NAUCHNIE

LABS QUALITY PRODUCTS

SRB-190 biological microscope with build in rechargeable batteries



Description

		SRB-190M	SRB-190B	SRB-190T
Optical System	Finite Optical System	•	•	•
Viewing Head	Monocular Head, Inclined at 30° , 360°Rotatable	•		
	Seidentopf Binocular Viewing Head, Inclined at 30° , 360°Rotatable		•	
	Seidentopf Trinocular Head, Inclined at 30° , 360°Rotatable			•
Eyepiece	Wide Field Eyepiece WF10X/18mm	•	•	•
	Wide Field Eyepiece WF16X/13mm	○	○	○
	Wide Field Eyepiece WF20X/11mm	○	○	○
Objective	4X 10X 40Xs 100Xs(oil) Finite 195 Series Achromatic Objectives 4X 10X 40Xs 100Xs(oil)	•	•	•
	Finite Plan Achromatic CF Objectives 4X 10X 40Xs 100Xs(oil)	○	○	○
Nosepiece	Quadruple Nosepiece	•	•	•
Condenser	Abbe Condenser NA1.25 with Diaphragm, Filter	•	•	•
Focusing System	Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Coarse Stroke36mm	•	•	•
Stage	Double Layer Mechanical Stage 110x125mm	•	•	•
Illumination	3W/LED Adjustment, rechargeable battery	•	•	•