

BIO BASIC INC.

160 TORBAY ROAD, MARKHAM ONTARIO L3R 1G6 CANADA
TEL: (905) 474 4493, (800) 313 7224 FAX: (905) 474 5794

MATERIAL SAFETY DATA SHEET

REQUIRED UNDER SAFETY AND HEALTH REGULATION FOR SHIP REPAIRING

MANUFACTURER'S NAME: BIO BASIC INC.
TEL: (905) 474 4493
FAX: (905) 474 5794
FOR CHEMICAL EMERGENCY: (416) 995 9730
ADDRESS: 160 TORBAY ROAD, MARKHAM ONTARIO L3R 1G6 CANADA

DATE PREPARED: JULY 25, 2007

SECTION 1 - PRODUCT AND COMPANY INFORMATION

PRODUCT NAME ANTHRAQUINONE
PRODUCT NUMBER AB0087

SECTION 2 - COMPOSITION/INFORMATION ON INGREDIENT

SUBSTANCE NAME	CAS #	SARA 313
ANTHRAQUINONE	84-65-1	NO

FORMULA C₁₄H₈O₂

SYNONYMS ANTHRACENE, 9,10-DIHYDRO-9,10-DIOXO- *
9,10-ANTHRACENEDIONE * 9,10-ANTHRACHINON (CZECH)
* ANTHRADIONE * 9,10-ANTHRAQUINONE * CORBIT *
9,10-DIOXOANTHRACENE * HOELITE * MORKIT
RTECS NUMBER: CB4725000

SECTION 3 - HAZARDS IDENTIFICATION

EMERGENCY OVERVIEW

HARMFUL.

IRRITATING TO EYES, RESPIRATORY SYSTEM AND SKIN. LIMITED EVIDENCE OF A CARCINOGENIC EFFECT. MAY CAUSE SENSITIZATION BY SKIN CONTACT. POSSIBLE CARCINOGEN (US). TARGET ORGAN(S): KIDNEYS. LIVER. BLADDER.

HMIS RATING

HEALTH: 3*

FLAMMABILITY: 1

REACTIVITY: 0

NFPA RATING

HEALTH: 3

FLAMMABILITY: 1

REACTIVITY: 0

*ADDITIONAL CHRONIC HAZARDS PRESENT.

FOR ADDITIONAL INFORMATION ON TOXICITY, PLEASE REFER TO SECTION 11.

SECTION 4 - FIRST AID MEASURES

ORAL EXPOSURE

IF SWALLOWED, WASH OUT MOUTH WITH WATER PROVIDED PERSON IS CONSCIOUS. CALL A PHYSICIAN.

INHALATION EXPOSURE

IF INHALED, REMOVE TO FRESH AIR. IF BREATHING BECOMES DIFFICULT, CALL A PHYSICIAN.

DERMAL EXPOSURE

IN CASE OF SKIN CONTACT, FLUSH WITH COPIOUS AMOUNTS OF WATER FOR AT LEAST 15 MINUTES. REMOVE CONTAMINATED CLOTHING AND SHOES. CALL A PHYSICIAN.

EYE EXPOSURE

IN CASE OF CONTACT WITH EYES, FLUSH WITH COPIOUS AMOUNTS OF WATER FOR AT LEAST 15 MINUTES. ASSURE ADEQUATE FLUSHING BY SEPARATING THE EYELIDS WITH FINGERS. CALL A PHYSICIAN.

SECTION 5 - FIRE FIGHTING MEASURES

FLASH POINT

365 °F 185 °C METHOD: CLOSED CUP

AUTOIGNITION TEMP

N/A

FLAMMABILITY

N/A

EXTINGUISHING MEDIA

SUITABLE: WATER SPRAY. CARBON DIOXIDE, DRY CHEMICAL POWDER, OR APPROPRIATE FOAM.

FIRE FIGHTING

PROTECTIVE EQUIPMENT: WEAR SELF-CONTAINED BREATHING APPARATUS AND PROTECTIVE CLOTHING TO PREVENT CONTACT WITH SKIN AND EYES.

SPECIFIC HAZARD(S): EMITS TOXIC FUMES UNDER FIRE CONDITIONS.

SECTION 6 - ACCIDENTAL RELEASE MEASURES

PROCEDURE(S) OF PERSONAL PRECAUTION(S)

WEAR RESPIRATOR, CHEMICAL SAFETY GOGGLES, RUBBER BOOTS, AND HEAVY RUBBER GLOVES.

METHODS FOR CLEANING UP

SWEEP UP, PLACE IN A BAG AND HOLD FOR WASTE DISPOSAL. AVOID RAISING DUST. VENTILATE AREA AND WASH SPILL SITE AFTER MATERIAL PICKUP IS COMPLETE.

SECTION 7 - HANDLING AND STORAGE

HANDLING

USER EXPOSURE: AVOID BREATHING DUST. AVOID CONTACT WITH EYES, SKIN, AND CLOTHING. AVOID PROLONGED OR REPEATED EXPOSURE.

STORAGE

SUITABLE: KEEP TIGHTLY CLOSED.

SECTION 8 - EXPOSURE CONTROLS / PPE

ENGINEERING CONTROLS

MECHANICAL EXHAUST REQUIRED.

PERSONAL PROTECTIVE EQUIPMENT

RESPIRATORY: USE RESPIRATORS AND COMPONENTS TESTED AND APPROVED UNDER APPROPRIATE GOVERNMENT STANDARDS SUCH AS NIOSH (US) OR CEN (EU). WHERE RISK ASSESSMENT SHOWS AIR-PURIFYING RESPIRATORS ARE APPROPRIATE USE A DUST MASK TYPE N95 (US) OR TYPE P1 (EN 143) RESPIRATOR.

OTHER: WEAR APPROPRIATE GOVERNMENT APPROVED RESPIRATOR, CHEMICAL-RESISTANT GLOVES, SAFETY GOGGLES, OTHER PROTECTIVE CLOTHING.

GENERAL HYGIENE MEASURES

WASH THOROUGHLY AFTER HANDLING. WASH CONTAMINATED CLOTHING BEFORE REUSE.

SECTION 9 - PHYSICAL/CHEMICAL PROPERTIES

APPEARANCE COLOR: LIGHT YELLOW

FORM: POWDER

PROPERTY VALUE AT TEMPERATURE OR PRESSURE

MOLECULAR WEIGHT 208.22 AMU

PH N/A

BP/BP RANGE 379.0 - 381.0 °C 760 MMHG
MP/MP RANGE 283 °C
FREEZING POINT N/A
VAPOR PRESSURE 1 MMHG 190 °C
VAPOR DENSITY 7.16 G/L
SATURATED VAPOR CONC. N/A
SG/DENSITY N/A
BULK DENSITY N/A
ODOR THRESHOLD N/A
VOLATILE% N/A
VOC CONTENT N/A
WATER CONTENT N/A
SOLVENT CONTENT N/A
EVAPORATION RATE N/A
VISCOSITY N/A
SURFACE TENSION N/A
PARTITION COEFFICIENT N/A
DECOMPOSITION TEMP. N/A
FLASH POINT 365 °F 185 °C METHOD: CLOSED CUP
EXPLOSION LIMITS N/A
FLAMMABILITY N/A
AUTOIGNITION TEMP N/A
REFRACTIVE INDEX N/A
OPTICAL ROTATION N/A
MISCELLANEOUS DATA N/A
SOLUBILITY N/A
N/A = NOT AVAILABLE

SECTION 10 - STABILITY AND REACTIVITY

STABILITY

STABLE: STABLE.

MATERIALS TO AVOID: OXIDIZING AGENTS.

HAZARDOUS DECOMPOSITION PRODUCTS

HAZARDOUS DECOMPOSITION PRODUCTS: CARBON MONOXIDE, CARBON DIOXIDE.

HAZARDOUS POLYMERIZATION

HAZARDOUS POLYMERIZATION: WILL NOT OCCUR

SECTION 11 - TOXICOLOGICAL INFORMATION

ROUTE OF EXPOSURE

SKIN CONTACT: CAUSES SKIN IRRITATION.

SKIN ABSORPTION: MAY BE HARMFUL IF ABSORBED THROUGH THE SKIN.

EYE CONTACT: CAUSES EYE IRRITATION.

INHALATION: MATERIAL IS IRRITATING TO MUCOUS MEMBRANES AND UPPER

RESPIRATORY TRACT. MAY BE HARMFUL IF INHALED.

INGESTION: MAY BE HARMFUL IF SWALLOWED.

SENSITIZATION

SKIN: MAY CAUSE ALLERGIC SKIN REACTION.

TARGET ORGAN(S) OR SYSTEM(S)

BLADDER. LIVER. KIDNEYS.

SIGNS AND SYMPTOMS OF EXPOSURE

TO THE BEST OF OUR KNOWLEDGE, THE CHEMICAL, PHYSICAL, AND

TOXICOLOGICAL PROPERTIES HAVE NOT BEEN THOROUGHLY INVESTIGATED.

TOXICITY DATA

INHALATION

RAT

> 1,300 MG/M3

LC50

SKIN

RAT

> 1000 MG/KG
LD50
INTRAPERITONEAL
RAT

3500 MG/KG
LD50
ORAL
MOUSE

> 5000 MG/KG
LD50

CHRONIC EXPOSURE - CARCINOGEN

RESULT: THIS PRODUCT IS OR CONTAINS A COMPONENT THAT HAS BEEN REPORTED TO BE POSSIBLY CARCINOGENIC BASED ON ITS IARC, ACGIH, NTP, OR EPA CLASSIFICATION.

NTP CARCINOGEN LIST

RATING: SOME EVIDENCE.

SPECIES: RAT, MALE

ROUTE: FEED

CHRONIC EXPOSURE - MUTAGEN

SPECIES: MOUSE

ROUTE: FEED

MUTATION TEST: MICRONUCLEUS TEST

RESULT: LABORATORY EXPERIMENTS HAVE SHOWN MUTAGENIC EFFECTS.

SPECIES: MOUSE

ROUTE: FEED

CELL TYPE: BONE MARROW

MUTATION TEST: MICRONUCLEUS TEST

SPECIES: MOUSE

ROUTE: INTRAPERITONEAL

DOSE: 250 MG/KG

MUTATION TEST: DNA DAMAGE

SPECIES: HAMSTER

DOSE: 12500 UG/L

CELL TYPE: EMBRYO

MUTATION TEST: MICRONUCLEUS TEST

SECTION 12 - ECOLOGICAL INFORMATION

NO DATA AVAILABLE.

SECTION 13 - DISPOSAL CONSIDERATIONS

APPROPRIATE METHOD OF DISPOSAL OF SUBSTANCE OR PREPARATION CONTACT A LICENSED PROFESSIONAL WASTE DISPOSAL SERVICE TO DISPOSE OF THIS MATERIAL. DISSOLVE OR MIX THE MATERIAL WITH A COMBUSTIBLE SOLVENT AND BURN IN A CHEMICAL INCINERATOR EQUIPPED WITH AN AFTERBURNER AND SCRUBBER. OBSERVE ALL FEDERAL, STATE, AND LOCAL ENVIRONMENTAL REGULATIONS.

SECTION 14 - TRANSPORT INFORMATION

DOT

PROPER SHIPPING NAME: NONE

NON-HAZARDOUS FOR TRANSPORT: THIS SUBSTANCE IS CONSIDERED TO BE NON-HAZARDOUS FOR TRANSPORT.

IATA

NON-HAZARDOUS FOR AIR TRANSPORT: NON-HAZARDOUS FOR AIR TRANSPORT.

SECTION 15 - REGULATORY INFORMATION

EU ADDITIONAL CLASSIFICATION

SYMBOL OF DANGER: XI

INDICATION OF DANGER: IRRITANT.

R: 43

RISK STATEMENTS: MAY CAUSE SENSITIZATION BY SKIN CONTACT.

S: 36/37

SAFETY STATEMENTS: WEAR SUITABLE PROTECTIVE CLOTHING AND GLOVES.

US CLASSIFICATION AND LABEL TEXT

INDICATION OF DANGER: HARMFUL.

RISK STATEMENTS: IRRITATING TO EYES, RESPIRATORY SYSTEM AND SKIN. LIMITED EVIDENCE OF A CARCINOGENIC EFFECT. MAY CAUSE SENSITIZATION BY SKIN CONTACT.

SAFETY STATEMENTS: IN CASE OF CONTACT WITH EYES, RINSE IMMEDIATELY WITH PLENTY OF WATER AND SEEK MEDICAL ADVICE. WEAR SUITABLE PROTECTIVE CLOTHING, GLOVES, AND EYE/FACE PROTECTION. IN CASE OF ACCIDENT OR IF YOU FEEL UNWELL, SEEK MEDICAL ADVICE IMMEDIATELY (SHOW THE LABEL WHERE POSSIBLE).

US STATEMENTS: POSSIBLE CARCINOGEN (US). TARGET ORGAN(S): KIDNEYS. LIVER. BLADDER.

UNITED STATES REGULATORY INFORMATION

SARA LISTED: NO

TSCA INVENTORY ITEM: YES

CANADA REGULATORY INFORMATION

WHMIS CLASSIFICATION: THIS PRODUCT HAS BEEN CLASSIFIED IN ACCORDANCE WITH THE HAZARD CRITERIA OF THE CPR, AND THE MSDS CONTAINS ALL THE INFORMATION REQUIRED BY THE CPR.

DSL: YES

NDSL: NO

SECTION 16 - OTHER INFORMATION

DISCLAIMER

FOR R&D USE ONLY. NOT FOR DRUG, HOUSEHOLD OR OTHER USES.

WARRANTY

THE ABOVE INFORMATION IS BELIEVED TO BE CORRECT BUT DOES NOT PURPORT TO BE ALL INCLUSIVE AND SHALL BE USED ONLY AS A GUIDE. THE INFORMATION IN THIS DOCUMENT IS BASED ON THE PRESENT STATE OF OUR KNOWLEDGE AND IS APPLICABLE TO THE PRODUCT WITH REGARD TO APPROPRIATE SAFETY PRECAUTIONS. IT DOES NOT REPRESENT ANY GUARANTEE OF THE PROPERTIES OF THE PRODUCT. BIOBASIC, SHALL NOT BE HELD LIABLE FOR ANY DAMAGE RESULTING FROM HANDLING OR FROM CONTACT WITH THE ABOVE PRODUCT. SEE REVERSE SIDE OF INVOICE OR PACKING SLIP FOR ADDITIONAL TERMS AND CONDITIONS OF SALE. COPYRIGHT 2007 BIOBASIC. LICENSE GRANTED TO MAKE UNLIMITED PAPER COPIES FOR INTERNAL USE ONLY.