



# SAFETY DATA SHEET

[Required under safety and health regulations for shipping and handling]

Version: 2020  
Date Updated: October 28, 2020

## SECTION 1. ----- PRODUCT AND COMPANY IDENTIFICATION -----

**Product Name** D-Arabinose  
**Product Code(s)** AB0071D  
**Recommended Use** For Laboratory Research Use Only  
Not for Human or Animal Drug Use

**Supplier** Bio Basic Inc.  
**Address** 20 Konrad Crescent, Markham, Ontario,  
Canada, L3R 8T4  
**Telephone** (905) 474 4493  
**Fax** (905) 474 5794  
**For Chemical Emergency Phone#** (416) 995 9730

## SECTION 2. ----- HAZARDS IDENTIFICATION -----

### **Classification of the substance or mixture**

Not a hazardous substance or mixture.

### **GHS Label elements, including precautionary statements**

Not a hazardous substance or mixture

**Hazards not otherwise classified (HNOC) or not covered by GHS** - none

### **Potential Health Effects**

**Inhalation** May be harmful if inhaled. May cause respiratory tract irritation.  
**Skin** May be harmful if absorbed through skin. May cause skin irritation.  
**Eyes** May cause eye irritation.  
**Ingestion** May be harmful if swallowed.

## SECTION 3. ----- COMPOSITION/INFORMATION ON INGREDIENTS -----

Chemical Name	EC No.	CAS-No	Weight %
D-Arabinose	233-708-5	10323-20-3	<100

## SECTION 4. ----- FIRST-AID MEASURES -----

### **If inhaled**

If breathed in, move person into fresh air. If not breathing, give artificial respiration.

### **In case of skin contact**

Wash off with soap and plenty of water.

### **In case of eye contact**

Flush eyes with water as a precaution.

### **If swallowed**

Never give anything by mouth to an unconscious person. Rinse mouth with water.

**Most important symptoms and effects, both acute and delayed**

The most important known symptoms and effects are described in the labelling (see section 2) and/or in section 11

**Indication of any immediate medical attention and special treatment needed**

No data available

**SECTION 5. - - - - - FIRE FIGHTING MEASURES - - - - -**

**Suitable extinguishing media**

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

**Special hazards arising from the substance or mixture**

Carbon oxides

**Special protective equipment for firefighters**

Wear self contained breathing apparatus for fire fighting if necessary.

**Hazardous combustion products**

**Explosion data - sensitivity to mechanical impact**

no data available

**Explosion data - sensitivity to static discharge**

no data available

**SECTION 6. - - - - - ACCIDENTAL RELEASE MEASURES- - - - -**

**Personal precautions, protective equipment and emergency procedures**

Avoid dust formation. Avoid breathing vapours, mist or gas. For personal protection see section 8.

**Environmental precautions**

Do not let product enter drains.

**Methods and materials for containment and cleaning up**

Sweep up and shovel. Keep in suitable, closed containers for disposal.

**Reference to other sections**

For disposal see section 13.

**SECTION 7. - - - - - HANDLING AND STORAGE- - - - -**

**Precautions for safe handling**

Provide appropriate exhaust ventilation at places where dust is formed. For precautions see section 2.

**Conditions for safe storage**

Keep container tightly closed in a dry and well-ventilated place.

Keep in a dry place.

Storage class (TRGS 510): 13: Non Combustible Solids

**Specific end use(s)**

Apart from the uses mentioned in section 1 no other specific uses are stipulated

**SECTION 8. - - - - - EXPOSURE CONTROLS/PERSONAL PROTECTION- - - - -**

**Appropriate engineering controls**

General industrial hygiene practice.

### Personal protective equipment

#### Respiratory protection

Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type N95 (US) or type P1 (EN 143) dust masks. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

#### Skin and Hand protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Full contact

Material: Nitrile rubber

Minimum layer thickness: 0.11 mm

Break through time: 480 min

Splash contact Material: Nitrile rubber

Minimum layer thickness: 0.11 mm

Break through time: 480 min

If used in solution, or mixed with other substances, and under conditions which differ from EN 374, contact the supplier of the CE approved gloves. This recommendation is advisory only and must be evaluated by an industrial hygienist and safety officer familiar with the specific situation of anticipated use by our customers. It should not be construed as offering an approval for any specific use scenario.

#### Eye protection

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

#### Body protection

Choose body protection in relation to its type, to the concentration and amount of dangerous substances, and to the specific work-place., The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

#### Control of environmental exposure

Do not let product enter drains.

## SECTION 9. - - - - - PHYSICAL AND CHEMICAL PROPERTIES - - - - -

### Appearance

Form	powder
Colour	white

### Safety data

pH	no data available
Melting point/freezing point	Melting point/range: 162 - 164 °C (324 - 327 °F) - lit.
Boiling point	no data available
Flash point	no data available
Ignition temperature	no data available

Auto-ignition temperature	no data available
Lower explosion limit	no data available
Upper explosion limit	no data available
Vapour pressure	no data available
Density	no data available
Water solubility	No data available
Partition coefficient: n-octanol/water	no data available
Relative vapour density	No data available
Odour	no data available
Odour Threshold	no data available
Evaporation rate	no data available

## SECTION 10. -----STABILITY AND REACTIVITY -----

### **Reactivity**

No data available

### **Chemical stability**

Stable under recommended storage conditions.

### **Possibility of hazardous reactions**

no data available

### **Conditions to avoid**

No data available

### **Incompatible Materials**

Strong oxidizing agents

### **Hazardous decomposition products**

Hazardous decomposition products formed under fire conditions. - Carbon oxides

Other decomposition products - No data available

In the event of fire: see section 5

## SECTION 11. ----- TOXICOLOGICAL INFORMATION -----

### **Acute toxicity**

#### **Oral LD50**

No data available

#### **Inhalation LC50**

no data available

#### **Dermal LD50**

no data available

### **Other information on acute toxicity**

no data available

**Skin corrosion/irritation**

No data available

**Serious eye damage/eye irritation**

no data available

**Respiratory or skin sensitisation**

no data available

**Germ cell mutagenicity**

no data available

**Carcinogenicity**

This product is or contains a component that has been reported to be possibly carcinogenic based on its IARC, ACGIH, NTP, or EPA classification.

Limited evidence of carcinogenicity in animal studies

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

**Reproductive toxicity**

no data available

**Teratogenicity**

no data available

**Specific target organ toxicity - single exposure (Globally Harmonized System)**

no data available

**Specific target organ toxicity - repeated exposure (Globally Harmonized System)**

no data available

**Aspiration hazard**

no data available

**Synergistic effects**

no data available

**Additional Information**

RTECS: Not available

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

**SECTION 12. ----- ECOLOGICAL INFORMATION -----**

**Toxicity**

no data available

**Persistence and degradability**

no data available

**Bioaccumulative potential**

no data available

**Mobility in soil**

no data available

**PBT and vPvB assessment**

PBT/vPvB assessment not available as chemical safety assessment not required/not conducted

**Other adverse effects**

no data available

**SECTION 13. ----- DISPOSAL CONSIDERATIONS -----****Product**

Offer surplus and non-recyclable solutions to a licensed disposal company.

**Contaminated packaging**

Dispose of as unused product.

**SECTION 14. ----- TRANSPORT INFORMATION -----****DOT (US)**

Not dangerous goods

**IMDG**

Not dangerous goods

**IATA**

Not dangerous goods

**SECTION 15. ----- REGULATORY INFORMATION -----**

This product has been classified in accordance with the hazard criteria of the Hazardous Products Regulations (HPR) and the SDS contains all the information required by the HPR.

**SECTION 16. ----- OTHER INFORMATION -----**

Further information: no limited for paper copy, just for internal uses.

For research use only. Not intended for human or animal diagnostic or therapeutic uses.

**Disclaimer**

The information provided on this MSDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

**Issuing Date:** 28-Oct-2020

**End of SDS**