



SAFETY DATA SHEET

[Required under safety and health regulations for shipping and handling]

Version: 2020

Date Updated: October 28, 2020

SECTION 1. ----- PRODUCT AND COMPANY IDENTIFICATION-----

Product Name Calcium chloride, dihydrate
 Product Code(s) CD0050
 Recommended Use For Laboratory Research Use Only
 Not for Human or Animal Drug Use


Supplier Bio Basic Inc.
 Address 20 Konrad Crescent, Markham, Ontario,
 Canada, L3R 8T4
 Telephone (905) 474 4493
 Fax (905) 474 5794
 For Chemical Emergency Phone# (416) 995 9730

SECTION 2. ----- HAZARDS IDENTIFICATION -----

GHS Classification

Eye irritation (Category 2A), H319
 For the full text of the H-Statements mentioned in this Section, see Section 16.

GHS Label elements, including precautionary statements

Pictogram 

Signal word Warning

Hazard statement(s)
 H319 Causes serious eye irritation.

Precautionary statement(s)
 P264 Wash skin thoroughly after handling.
 P280 Wear eye protection/ face protection.
 P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
 P337 + P313 If eye irritation persists: Get medical advice/ attention. Use personal protective equipment as required.

Hazards not otherwise classified (HNOC) or not covered by GHS - none

SECTION 3. ----- COMPOSITION/INFORMATION ON INGREDIENTS -----

Chemical Name	EC No.	CAS-No	Weight %
Calcium chloride, dihydrate		10035-04-8	<100

SECTION 4. ----- FIRST-AID MEASURES-----

General advice

Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

In case of skin contact

Wash off with soap and plenty of water. Consult a physician.

In case of eye contact

Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (see section 2) and/or in section 11

Indication of any immediate medical attention and special treatment needed

no data available

SECTION 5. ----- FIRE FIGHTING MEASURES -----**Suitable extinguishing media**

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special hazards arising from the substance or mixture

Hydrogen chloride gas, Calcium oxide
Nature of decomposition products not known

Special protective equipment for firefighters

Wear self-contained breathing apparatus for firefighting if necessary

Hazardous combustion products**Explosion data - sensitivity to mechanical impact**

no data available

Explosion data - sensitivity to static discharge

no data available

SECTION 6. ----- ACCIDENTAL RELEASE MEASURES-----**Personal precautions, protective equipment and emergency procedures**

Use personal protective equipment. Avoid dust formation. Avoid breathing vapours, mist or gas. Ensure adequate ventilation. Avoid breathing dust. For personal protection see section 8.

Environmental precautions

Do not let product enter drains.

Methods and materials for containment and cleaning up

Pick up and arrange disposal without creating dust. Sweep up and shovel. Keep in suitable, closed containers for disposal.
For disposal see section 13.

SECTION 7. ----- HANDLING AND STORAGE-----**Precautions for safe handling**

Avoid contact with skin and eyes. Avoid formation of dust and aerosols. Provide appropriate exhaust ventilation at places where dust is formed. For precautions see section 2.

Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place. Storage class (TRGS 510): 13: Non Combustible Solids

Specific end use(s)

Apart from the uses mentioned in section 1 no other specific uses are stipulated

SECTION 8. - - - - EXPOSURE CONTROLS/PERSONAL PROTECTION- - - -

Control parameters

Components with workplace control parameters

<i>Compon ents</i>	<i>CAS-No.</i>	<i>Value</i>	<i>Control parameters</i>	<i>Basis</i>
Calcium chloride dihydrate	10035-04- 8	TWA	5 mg/m3	Ontario Table of Occupational Exposure Limits made under the Occupational Health and Safety Act.

Personal protective equipment

Specific engineering controls

Handle in accordance with good industrial hygiene and safety practice. Wash hands before breaks and at the end of workday.

Respiratory protection

For nuisance exposures use type P95 (US) or type P1 (EU EN 143) particle respirator. For higher level protection use type OV/AG/P99 (US) or type ABEK-P2 (EU EN 143) respirator cartridges. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Full contact

Material: Nitrile rubber

Minimum layer thickness: 0.11 mm

Break through time: 480 min

Material tested: Dermatril® (KCL 740 / Aldrich Z677272, Size M)

Splash contact

Material: Nitrile rubber

Minimum layer thickness: 0.11 mm

Break through time: 480 min

Material tested: Dermatril® (KCL 740 / Aldrich Z677272, Size M)

data source: KCL GmbH, D-36124 Eichenzell, phone +49 (0)6659 87300, e-mail sales@kcl.de, test method: EN374

If used in solution, or mixed with other substances, and under conditions which differ from EN 374, contact the supplier of the CE approved gloves. This recommendation is advisory only and must be evaluated by an industrial hygienist and safety officer familiar with the specific situation of anticipated use by our customers. It should not be construed as offering an approval for any specific use scenario.

Eye protection

Safety glasses with side-shields conforming to EN166 Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Skin and body protection

Impervious clothing, The type of protective equipment must be selected according to the concentration

and amount of the dangerous substance at the specific workplace.

SECTION 9. ----- PHYSICAL AND CHEMICAL PROPERTIES -----

Appearance

Form	crystalline
Colour	white

Safety data

pH	5.0 - 8 at 147 g/l at 25 °C (77 °F)
Melting point/freezing point	Melting point/range: 176 °C (349 °F) - dec.
Boiling point	No data available
Flash point	no data available
Ignition temperature	no data available
Auto-ignition temperature	no data available
Lower explosion limit	no data available
Upper explosion limit	no data available
Vapour pressure	0.01 hPa at 20 °C (68 °F)
Relative Density	1.850 g/cm ³
Water solubility	147 g/l at 20 °C (68 °F) - completely soluble
Partition coefficient: n-octanol/water	no data available
Relative vapour density	no data available
Odour	no data available
Odour Threshold	no data available
Evaporation rate	no data available

SECTION 10. ----- STABILITY AND REACTIVITY -----

Reactivity

No data available

Chemical stability

Stable under recommended storage conditions.

Possibility of hazardous reactions

no data available

Conditions to avoid

No data available

Incompatible materials

Strong acids, Borane/boron oxides, Zinc, Calcium oxide, Methyl vinyl ether, Calcium chloride is attacked by bromine trifluoride

Hazardous decomposition products

Hazardous decomposition products formed under fire conditions. - Hydrogen chloride gas,
Calcium oxide
Hazardous decomposition products formed under fire conditions. - Nature of decomposition
products not known.
Other decomposition products - No data available In the event of fire: see section 5

SECTION 11. ----- TOXICOLOGICAL INFORMATION -----

Acute toxicity

Oral LD50

LD50 Oral - Rabbit - male - 500 - 1,000 mg/kg (OECD Test Guideline 401) Remarks: (anhydrous substance)

Inhalation LC50

Irritating to respiratory system.

Dermal LD50

LD50 Dermal - Rabbit - male and female - > 5,000 mg/kg
Remarks: (anhydrous substance) (ECHA)

Other information on acute toxicity

no data available

Skin corrosion/irritation

Skin - Rabbit

Result: No skin irritation - 4 h
(OECD Test Guideline 404)

Remarks: (anhydrous substance)

Serious eye damage/eye irritation

Eyes - Rabbit

Result: Moderate eye irritation
(OECD Test Guideline 405)

Remarks: (anhydrous substance)

Respiratory or skin sensitisation

no data available

Germ cell mutagenicity

Ames test

S. typhimurium

Result: negative

Mutagenicity (mammal cell test): chromosome aberration.

Chinese hamster fibroblasts

Result: negative

Carcinogenicity

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

Reproductive toxicity

no data available

Teratogenicity

no data available

Specific target organ toxicity - single exposure (Globally Harmonized System)

Acute oral toxicity - After uptake of large quantities: Stomach/intestinal disorders, Nausea

Specific target organ toxicity - repeated exposure (Globally Harmonized System)

no data available

Aspiration hazard

no data available

Synergistic effects

no data available

Additional Information

RTECS: EV9810000

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

SECTION 12. ----- ECOLOGICAL INFORMATION -----

Toxicity

Toxicity to fish static test LC50 - Pimephales promelas (fathead minnow) - 4,630 mg/l - 96 h
(US-EPA)
Remarks: (anhydrous substance)

Toxicity to daphnia and other aquatic invertebrates static test

EC50 - Daphnia magna (Water flea) - 2,400 mg/l - 48 h
(OECD Test Guideline 202)
Remarks: (anhydrous substance)

Toxicity to algae EC50 - Pseudokirchneriella subcapitata - 2,900 mg/l - 72 h
(OECD Test Guideline 201)
Remarks: (anhydrous substance)

Persistence and degradability

The methods for determining biodegradability are not applicable to inorganic substances.

Bioaccumulative potential

no data available

Mobility in soil

no data available

PBT and vPvB assessment

PBT/vPvB assessment not available as chemical safety assessment not required/not conducted

Other adverse effects

no data available

SECTION 13. ----- DISPOSAL CONSIDERATIONS -----

Product

Offer surplus and non-recyclable solutions to a licensed disposal company. Contact a licensed professional waste disposal service to dispose of this material.

Contaminated packaging

Dispose of as unused product.

SECTION 14. ----- TRANSPORT INFORMATION -----

DOT (US)

Not dangerous goods

IMDG

Not dangerous goods

IATA

Not dangerous goods

SECTION 15. ----- REGULATORY INFORMATION -----

This product has been classified in accordance with the hazard criteria of the Hazardous Products Regulations (HPR) and the SDS contains all the information required by the HPR.

SECTION 16. ----- OTHER INFORMATION -----

Further information: no limited for paper copy, just for internal uses.
For research use only. Not intended for human or animal diagnostic or therapeutic uses.

Disclaimer

The information provided on this MSDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

Issuing Date: 28-Oct-2020

End of SDS