

## BIO BASIC INC.

160 TORBAY ROAD MARKHAM ONTARIO L3R 1G6 CANADA  
TEL: (905) 474 4493, (800) 313 7224 FAX: (905) 474 5794

### MATERIAL SAFETY DATA SHEET

REQUIRED UNDER SAFETY AND HEALTH REGULATION FOR SHIP REPAIRING  
MANUFACTURER'S NAME: BIO BASIC INC.

TEL: (905) 474 4493

FAX: (905) 474 5794

FOR CHEMICAL EMERGENCY: (416) 995 9730

ADDRESS: 160 TORBAY ROAD MARKHAM ONTARIO L3R 1G6 CANADA

**DATE UPDATED: MAY 2007**

#### SECTION 1. ----- CHEMICAL IDENTIFICATION -----

CODE: AB0060

NAME: AMMONIUM CITRATE, TRIBASIC

#### Section 2 - Composition/Information on Ingredient

Substance Name

AMMONIUM CITRATE TRIBASIC

CAS # 3458-72-8

SARA NO

Please refer to the formulation sheet for exact formula composition.

#### Section 3 - Hazards Identification

HMIS Rating

HEALTH: 2\*

FLAMMABILITY: 0

REACTIVITY: 3

NFPA Rating

HEALTH: 2

FLAMMABILITY: 0

REACTIVITY: 3

For additional information on toxicity, please refer to Section 11.

Irritant. Irritating to eyes, respiratory system and skin.

Harmful if swallowed.

Toxic. Toxic by inhalation, in contact with skin and if swallowed.

\* Additional chronic hazards present.

#### Section 4 - First Aid Measures

##### ORAL EXPOSURE

If swallowed, wash out mouth with water provided person is conscious. Call a physician immediately. Do not induce vomiting.

##### INHALATION EXPOSURE

If inhaled, remove to fresh air. If not breathing give artificial respiration. If breathing is difficult, give oxygen. Call a physician.

##### DERMAL EXPOSURE

In case of contact, immediately wash skin with copious amounts of water for at least 15 minutes. Remove contaminated clothing and shoes. Call a physician.

##### EYE EXPOSURE

In case of contact with eyes, flush with copious amounts of water for at least 15 minutes. Assure adequate flushing by separating the eyelids with fingers. Call a physician.

#### Section 5 - Fire Fighting Measures

##### FLASH POINT

N/A

##### EXPLOSION LIMITS

N/A

##### AUTOIGNITION TEMP

N/A

##### FLAMMABILITY

N/A

##### EXTINGUISHING MEDIA

Suitable: Carbon dioxide, dry chemical powder, or appropriate foam.

##### FIREFIGHTING

Protective Equipment: Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.

Specific Hazard(s): Emits toxic fumes under fire conditions.

#### Section 6 - Accidental Release Measures

##### METHODS FOR CLEANING UP

Cover with dry lime or soda ash, pick up, keep in a closed container, and hold for waste disposal. Ventilate area and wash spill site after material pickup is complete.

#### Section 7 - Handling and Storage

##### Handling

User Exposure: Do not breathe vapor. Do not get into eyes, on skin, on clothing. Avoid prolonged or repeated exposure.

##### PROCEDURE TO BE FOLLOWED IN CASE OF LEAK OR SPILL

Evacuate area.

##### Storage

Suitable: Keep tightly closed.

##### PROCEDURE(S) OF PERSONAL PRECAUTION(S)

Wear self-contained breathing apparatus, rubber boots, and heavy rubber gloves.

#### Section 8 - Exposure Controls / PPE

##### ENGINEERING CONTROLS

Use only in a chemical fume hood. Safety shower and eye bath.

##### PERSONAL PROTECTIVE EQUIPMENT

Respiratory: Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type N95 (US) or type P1 (EN 143) dust masks.

Hand: Compatible chemical-resistant gloves.

Eye: Chemical safety goggles.

##### GENERAL HYGIENE MEASURES

Wash thoroughly after handling.

#### Section 9 - Physical/Chemical Properties

Property Value At Temperature or Pressure

Molecular Weight N/A

pH N/A  
BP/BP Range N/A  
MP/MP Range N/A  
Freezing Point N/A  
Vapor Pressure N/A  
Vapor Density N/A  
Saturated Vapor Conc. N/A  
SG/Density N/A  
Odor Threshold N/A  
Volatile % N/A  
VOC Content N/A  
Water Content N/A  
Solvent Content N/A  
Evaporation Rate N/A  
Viscosity N/A  
Surface Tension N/A  
Partition Coefficient N/A  
Decomposition Temp. N/A  
Appearance Physical State: Liquid  
Flash Point N/A  
Explosion Limits N/A  
Bulk Density N/A  
Flammability N/A  
Autoignition Temp. N/A  
Refractive Index N/A  
Optical Rotation N/A  
Miscellaneous Data N/A  
Solubility N/A

N/A = Not Available

#### Section 10 - Stability and Reactivity

##### STABILITY

Stable: Stable.

Materials to Avoid: Strong oxidizing agents.

##### HAZARDOUS DECOMPOSITION PRODUCTS

Hazardous Decomposition Products: Carbon monoxide, Carbon dioxide, Sodium/sodium oxides, Phosphorous oxides, Potassium oxides.

##### HAZARDOUS POLYMERIZATION

Hazardous Polymerization: Will not occur.

#### Section 11 - Toxicological Information

##### ROUTE OF EXPOSURE

Skin Contact: Causes skin irritation.

Skin Absorption: May be harmful if absorbed through the skin.

Eye Contact: Causes eye irritation.

Inhalation: Material is irritating to mucous membranes and upper respiratory tract. May be harmful if inhaled.

Ingestion: May be harmful if swallowed.

##### SIGNS AND SYMPTOMS OF EXPOSURE

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

#### Section 12 - Ecological Information

No data available.

Section 13 - Disposal Considerations

APPROPRIATE METHOD OF DISPOSAL OF SUBSTANCE OR PREPARATION

Contact a licensed professional waste disposal service to dispose of this material. Dissolve or mix the material with a combustible solvent and burn in a chemical incinerator equipped with an afterburner and scrubber. Observe all federal, state, and local environmental regulations.

Section 14 - Transport Information

DOT

Proper Shipping Name: None

Non-Hazardous for Transport: This substance is considered to be non-hazardous for transport.

IATA Non-Hazardous for Air Transport: Non-hazardous for air transport.

Section 15 - Regulatory Information

EU ADDITIONAL CLASSIFICATION

S: 24/25

Safety Statements: Avoid contact with skin and eyes.

US CLASSIFICATION AND LABEL TEXT

Indication of Danger: Irritant.

Risk Statements: Irritating to eyes, respiratory system and skin.

Safety Statements: Wear suitable protective clothing. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

UNITED STATES REGULATORY INFORMATION

SARA LISTED: No

CANADA REGULATORY INFORMATION

WHMIS Classification: This product has been classified in accordance with the hazard criteria of the CPR, and the MSDS contains all the information required by the CPR.

DSL: No

NDSL: No

Section 16 - Other Information

DISCLAIMER

For R&D use only. Not for drug, household or other uses.

WARRANTY

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of this product.

THE ABOVE INFORMATION IS BELIEVED TO BE CORRECT BUT DOES NOT PURPORT TO BE ALL INCLUSIVE AND SHALL BE USED ONLY AS A GUIDE. BIO BASIC INC. SHALL NOT BE HELD LIABLE FOR ANY DAMAGE RESULTING FROM HANDLING OR FROM CONTACT WITH THE ABOVE PRODUCT. SEE REVERSE SIDE OF INVOICE OR PACKING SLIP FOR ADDITIONAL TERMS AND CONDITIONS OF SALE.