

OVERVIEW

ZiXpress-32[®] (ZP01013)

Helps you to purify 1 to

32 samples

simultaneously per run



Specification

| | | ZiXpress 32 | ZiXpress 64 |
|-----------------------|-------------------|---------------------------|--------------------|
| GENERAL SPECIFICATION | Dimensions | 438 x 360 x 499 mm | 722 x 360 x 490 mm |
| | Weight | 25 kg | 38kg |
| | Throughput | 32 tests/run | 64 tests/run |
| | Power Supply | AC 100-240V, 400VA, UPS | |
| | Display | 7" inch color Touch Panel | |
| | Cat. No. | ZP01013 | ZP01014 |
| EXTRACTION FEATURES | Sample Volume | 200 µl | |
| | Elution Volume | 100-1000 µl | |
| | Processing Time | 30-45min | |
| TECHNOLOGY AND SAFETY | UV Decotamination | ✓ | ✓ |
| | Door Sensor | ✓ | ✓ |
| | Bottom Touch | ✓ | ✓ |