

# SAFETY DATA SHEET

[Required under safety and health regulations for shipping and handling]

Version: 2021

Date Updated: Apr 01, 2021

#### SECTION 1. - - - - - PRODUCT AND COMPANY IDENTIFICATION - - - - - -

Product Name Cycloheximide

Product Code(s) CB1351

Recommended Use For Laboratory Research Use Only

Not for Human or Animal Drug Use

**Supplier** Bio Basic Inc.

Address 20 Konrad Crescent, Markham, Ontario,

Canada, L3R 8T4

 Telephone
 (905) 474 4493

 Fax
 (905) 474 5794

 For Chemical Emergency Phone#
 (416) 995 9730

#### SECTION 2. ----- HAZARDS IDENTIFICATION -----

## GHS Classification in accordance with Hazardous Products Regulations (HPR) (SOR/2015-17)

Acute toxicity, Oral (Category 2), H300 Germ cell mutagenicity (Category 2), H341 Reproductive toxicity (Category 1B), H360

Short-term (acute) aquatic hazard (Category 2), H401 Long-term

(chronic) aquatic hazard (Category 2), H411

For the full text of the H-Statements mentioned in this Section, see Section 16.

# GHS Label elements, including precautionary statements

Pictogram

Signal word Danger

Hazard statement(s)

H300 Fatal if swallowed.

H341 Suspected of causing genetic defects.
H360 May damage fertility or the unborn child.
H411 Toxic to aquatic life with long lasting effects.

Precautionary statement(s)

P201 Obtain special instructions before use.

P202 Do not handle until all safety precautions have been read and

understood.

P264 Wash skin thoroughly after handling.

P270 Do not eat, drink or smoke when using this product.

P273 Avoid release to the environment.

P280 Wear protective gloves/ protective clothing/ eye protection/ face

protection.

P301 + P310 + P330 IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

Rinse mouth.

P308 + P313 IF exposed or concerned: Get medical advice/ attention.

P391 Collect spillage.

P405 Store locked up.

P501 Dispose of contents/ container to an approved waste disposal

plant.

# Hazards not otherwise classified (HNOC) or not covered by GHS

- none

# SECTION 3. ---- COMPOSITION/INFORMATION ON INGREDIENTS -----

Chemical Name	EC No.	CAS-No	Weight %
Cycloheximide	200-636-0	66-81-9	<100

For the full text of the H-Statements mentioned in this Section, see Section 16.

#### SECTION 4. - - - - - FIRST-AID MEASURES- - - - - -

#### General advice

Consult a physician. Show this material safety data sheet to the doctor in attendance. Move out of dangerous area.

#### If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

#### In case of skin contact

Wash off with soap and plenty of water. Take victim immediately to hospital. Consult a physician.

#### In case of eye contact

Flush eyes with water as a precaution.

#### If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

## Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (see section 2) and/or in section 11

#### Indication of any immediate medical attention and special treatment needed

No data available

## SECTION 5. ----- FIRE FIGHTING MEASURES -----

## Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

#### Special hazards arising from the substance or mixture

Carbon oxides, Nitrogen oxides (NOx)

# Advice for firefighters

Wear self-contained breathing apparatus for firefighting if necessary.

## **Further information**

No data available

## SECTION 6. ----- ACCIDENTAL RELEASE MEASURES-----

#### Personal precautions, protective equipment and emergency procedures

Wear respiratory protection. Avoid dust formation. Avoid breathing vapors, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas. Avoid breathing dust. For personal protection see section 8.

#### **Environmental precautions**

Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into the environment must be avoided.

### Methods and materials for containment and cleaning up

Pick up and arrange disposal without creating dust. Sweep up and shovel. Keep in suitable, closed containers for disposal.

#### Reference to other sections

For disposal see section 13.

## SECTION 7. ----- HANDLING AND STORAGE-----

#### Precautions for safe handling

Avoid contact with skin and eyes. Avoid formation of dust and aerosols. Provide appropriate exhaust ventilation at places where dust is formed. For precautions see section 2.

# Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place.

Keep in a dry place.

Storage class (TRGS 510): 6.1B: Non-combustible, acute toxic Cat. 1 and 2 / very toxic hazardous materials

## Specific end use(s)

Apart from the uses mentioned in section 1 no other specific uses are stipulated

#### SECTION 8. - - - - EXPOSURE CONTROLS/PERSONAL PROTECTION - - - -

## **Exposure controls**

## Appropriate engineering controls

Avoid contact with skin, eyes and clothing. Wash hands before breaks and immediately after handling the product.

# Personal protective equipment

# Eye/face protection

Face shield and safety glasses Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

## Skin protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Full contact

Material: Nitrile rubber

Minimum layer thickness: 0.11 mm Break

through timé: 480 min

Material tested:Dermatril® (KCL 740 / Aldrich Z677272, Size M)

Splash contact

Material: Nitrile rubber

Minimum layer thickness: 0.11 mm Break

through time: 480 min

Material tested: Dermatril® (KCL 740 / Aldrich Z677272, Size M)

data source: KCL GmbH, D-36124 Eichenzell, phone +49 (0)6659 87300, e-mail sales@kcl.de,

test method: EN374

If used in solution, or mixed with other substances, and under conditions which differ from EN 374, contact the supplier of the EC approved gloves. This recommendation is advisory only and

must be evaluated by an industrial hygienist and safety officer familiar with the specific situation of anticipated use by our customers. It should not be construed as offering an approval for any specific use scenario.

## **Body Protection**

Complete suit protecting against chemicals, The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

### Respiratory protection

Where risk assessment shows air-purifying respirators are appropriate use a full-face particle respirator type N100 (US) or type P3 (EN 143) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

## Control of environmental exposure

Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into the environment must be avoided.

## SECTION 9. ----- PHYSICAL AND CHEMICAL PROPERTIES -----

## Information on basic physical and chemical properties

iniormation on basic physical and chemical prop			
a)	Appearance	Form: powder Color: beige	
b)	Odor	No data available	
c)	Odor Threshold	No data available	
d)	рН	No data available	
e)	Melting point/freezing point	No data available	
f)	Initial boiling point and boiling range	No data available	
g)	Flash point	()No data available	
h)	Evaporation rate	No data available	
i)	Flammability (solid, gas)	No data available	
j)	Upper/lower flammability or explosive limits	No data available	
k) l)	Vapor pressure Vapor density	No data available No data available	
m)	Relative density	No data available	
n)	Water solubility	slightly soluble	
0)	Partition coefficient: log	g Pow:	
p)	Autoignition temperature	No data available	
q)	Decomposition	No data available	

temperature

Viscosity

QF26 Rev 2 4

No data available

s) Explosive properties No data available

t) Oxidizing properties No data available

## Other safety information

No data available

### SECTION 10. ------STABILITY AND REACTIVITY -----

## Reactivity

No data available

#### **Chemical stability**

Stable under recommended storage conditions.

## Possibility of hazardous reactions

No data available

#### Conditions to avoid

No data available

#### Incompatible materials

Bases, Strong oxidizing agents, Acid anhydrides, Acid chlorides

## **Hazardous decomposition products**

Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NOx) Other decomposition products - No data available In the event of fire: see section 5

# SECTION 11. ----- TOXICOLOGICAL INFORMATION -----

## **Acute toxicity**

No data available

Inhalation: No data available Dermal: No data available

#### Skin corrosion/irritation

No data available

# Serious eye damage/eye irritation

No data available

## Respiratory or skin sensitization

Prolonged or repeated exposure may cause allergic reactions in certain sensitive individuals.

#### Germ cell mutagenicity

Suspected of causing genetic defects.

## Carcinogenicity

IARC: No ingredient of this product present at levels greater than or equal to 0.1% is

identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No ingredient of this product present at levels greater than or equal to 0.1% is

identified as a carcinogen or potential carcinogen by ACGIH.

#### Reproductive toxicity

May damage the unborn child.

## Specific target organ toxicity - single exposure

No data available

#### Specific target organ toxicity - repeated exposure

No data available

## **Aspiration hazard**

No data available

#### **Additional Information**

RTECS: MA4375000

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

Liver - Irregularities - Based on Human Evidence

# SECTION 12. ----- ECOLOGICAL INFORMATION -----

### **Toxicity**

No data available

## Persistence and degradability

No data available

## Bioaccumulative potential

No data available

## Mobility in soil

No data available

#### Results of PBT and vPvB assessment

PBT/vPvB assessment not available as chemical safety assessment not required/not conducted

#### Other adverse effects

Toxic to aquatic life with long lasting effects.

An environmental hazard cannot be excluded in the event of unprofessional handling or disposal.

#### SECTION 13. ----- DISPOSAL CONSIDERATIONS -----

#### **Product**

Contact a licensed professional waste disposal service to dispose of this material. Dissolve or mix the material with a combustible solvent and burn in a chemical incinerator equipped with an afterburner and scrubber. Offer surplus and non-recyclable solutions to a licensed disposal company.

#### Contaminated packaging

Dispose of as unused product.

## SECTION 14. ----- TRANSPORT INFORMATION -----

# DOT (US)

UN number: 2811 Class: 6.1 Packing group: I Proper shipping name: Toxic solids, organic, n.o.s. (Cycloheximide)

Reportable Quantity (RQ): Poison Inhalation Hazard: No

#### **IMDG**

UN number: 2811 Class: 6.1 Packing group: I EMS-No: F-A, S-A Proper shipping name: TOXIC

SOLID, ORGANIC, N.O.S. (Cycloheximide)

Marine pollutant : yes

## **IATA**

UN number: 2811 Class: 6.1 Packing group: I Proper shipping name: Toxic solid, organic, n.o.s. (Cycloheximide)

## SECTION 15. ----- REGULATORY INFORMATION -----

This product has been classified in accordance with the hazard criteria of the Hazardous Products

Regulations (HPR) and the SDS contains all the information required by the HPR.

# SECTION 16. ----- OTHER INFORMATION-----

Further information: no limited for paper copy, just for internal uses. For research use only. Not intended for human or animal diagnostic or therapeutic uses.

## Disclaimer

The information provided on this MSDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

Issuing Date: 31-Mar-2021

**End of SDS**