

BIO BASIC INC.

160 TORBAY ROAD, MARKHAM ONTARIO L3R 1G6 CANADA
TEL: (905) 474 4493, (800) 313 7224 FAX: (905) 474 5794

MATERIAL SAFETY DATA SHEET

REQUIRED UNDER SAFETY AND HEALTH REGULATION FOR SHIP REPAIRING

MANUFACTURER'S NAME: BIO BASIC INC.
TEL: (905) 474 4493
FAX: (905) 474 5794
FOR CHEMICAL EMERGENCY: (416) 995 9730
ADDRESS: 160 TORBAY ROAD, MARKHAM ONTARIO L3R 1G6 CANADA

DATE PREPARED: JULY 30, 2007

SECTION 1 - PRODUCT AND COMPANY INFORMATION

PRODUCT NAME CREATININE
PRODUCT NUMBER CB0328

SECTION 2 - COMPOSITION/INFORMATION ON INGREDIENT

SUBSTANCE NAME	CAS #	SARA 313
CREATININE	60-27-5	NO
FORMULA C4H7N3O		

SECTION 3 - HAZARDS IDENTIFICATION

HMIS RATING
HEALTH: 0
FLAMMABILITY: 0
REACTIVITY: 0
NFPA RATING
HEALTH: 0
FLAMMABILITY: 0
REACTIVITY: 0

FOR ADDITIONAL INFORMATION ON TOXICITY, PLEASE REFER TO SECTION 11.

SECTION 4 - FIRST AID MEASURES

ORAL EXPOSURE

IF SWALLOWED, WASH OUT MOUTH WITH WATER PROVIDED PERSON IS CONSCIOUS. CALL A PHYSICIAN.

INHALATION EXPOSURE

IF INHALED, REMOVE TO FRESH AIR. IF BREATHING BECOMES DIFFICULT, CALL A PHYSICIAN.

DERMAL EXPOSURE

IN CASE OF CONTACT, IMMEDIATELY WASH SKIN WITH SOAP AND COPIOUS AMOUNTS OF WATER.

EYE EXPOSURE

IN CASE OF CONTACT WITH EYES, FLUSH WITH COPIOUS AMOUNTS OF WATER FOR AT LEAST 15 MINUTES. ASSURE ADEQUATE FLUSHING BY SEPARATING THE EYELIDS WITH FINGERS. CALL A PHYSICIAN.

SECTION 5 - FIRE FIGHTING MEASURES

FLASH POINT

554 °F 290 °C METHOD: CLOSED CUP

AUTOIGNITION TEMP

N/A

FLAMMABILITY

N/A

EXTINGUISHING MEDIA

SUITABLE: WATER SPRAY. CARBON DIOXIDE, DRY CHEMICAL POWDER, OR APPROPRIATE FOAM.

FIREFIGHTING

PROTECTIVE EQUIPMENT: WEAR SELF-CONTAINED BREATHING APPARATUS AND PROTECTIVE CLOTHING TO PREVENT CONTACT WITH SKIN AND EYES.

SPECIFIC HAZARD(S): EMITS TOXIC FUMES UNDER FIRE CONDITIONS.

SECTION 6 - ACCIDENTAL RELEASE MEASURES

PROCEDURE(S) OF PERSONAL PRECAUTION(S)

EXERCISE APPROPRIATE PRECAUTIONS TO MINIMIZE DIRECT CONTACT WITH SKIN OR EYES AND PREVENT INHALATION OF DUST.

METHODS FOR CLEANING UP

SWEEP UP, PLACE IN A BAG AND HOLD FOR WASTE DISPOSAL. AVOID RAISING DUST. VENTILATE AREA AND WASH SPILL SITE AFTER MATERIAL PICKUP IS COMPLETE.

SECTION 7 - HANDLING AND STORAGE

HANDLING

USER EXPOSURE: AVOID INHALATION. AVOID CONTACT WITH EYES, SKIN, AND CLOTHING. AVOID PROLONGED OR REPEATED EXPOSURE.

STORAGE

SUITABLE: KEEP TIGHTLY CLOSED.

SECTION 8 - EXPOSURE CONTROLS / PPE

ENGINEERING CONTROLS

SAFETY SHOWER AND EYE BATH. MECHANICAL EXHAUST REQUIRED.

PERSONAL PROTECTIVE EQUIPMENT

RESPIRATORY: USE RESPIRATORS AND COMPONENTS TESTED AND APPROVED UNDER APPROPRIATE GOVERNMENT STANDARDS SUCH AS NIOSH (US) OR CEN (EU). RESPIRATORY PROTECTION IS NOT REQUIRED. WHERE PROTECTION FROM NUISANCE LEVELS OF DUSTS ARE DESIRED, USE TYPE N95 (US) OR TYPE P1 (EN 143) DUST MASKS.

HAND: PROTECTIVE GLOVES.

EYE: CHEMICAL SAFETY GOGGLES.

GENERAL HYGIENE MEASURES

WASH THOROUGHLY AFTER HANDLING.

SECTION 9 - PHYSICAL/CHEMICAL PROPERTIES

APPEARANCE PHYSICAL STATE: SOLID

PROPERTY VALUE AT TEMPERATURE OR PRESSURE

MOLECULAR WEIGHT 113.12 AMU

PH 7.5 - 8.5 20 °C CONCENTRATION: 50

G/L

BP/BP RANGE N/A

MP/MP RANGE 255 °C

FREEZING POINT N/A

VAPOR PRESSURE N/A

VAPOR DENSITY N/A

SATURATED VAPOR CONC. N/A

SG/DENSITY N/A

BULK DENSITY N/A

ODOR THRESHOLD N/A

VOLATILE% N/A

VOC CONTENT N/A

WATER CONTENT N/A

SOLVENT CONTENT N/A

EVAPORATION RATE N/A

VISCOSITY N/A
SURFACE TENSION N/A
PARTITION COEFFICIENT N/A
DECOMPOSITION TEMP. N/A
FLASH POINT 554 °F 290 °C METHOD: CLOSED CUP
EXPLOSION LIMITS N/A
FLAMMABILITY N/A
AUTOIGNITION TEMP N/A
REFRACTIVE INDEX N/A
OPTICAL ROTATION N/A
MISCELLANEOUS DATA N/A
SOLUBILITY N/A
N/A = NOT AVAILABLE

SECTION 10 - STABILITY AND REACTIVITY

STABILITY

STABLE: STABLE.

MATERIALS TO AVOID: STRONG OXIDIZING AGENTS.

HAZARDOUS DECOMPOSITION PRODUCTS

HAZARDOUS DECOMPOSITION PRODUCTS: CARBON MONOXIDE, CARBON DIOXIDE,
NITROGEN OXIDES.

HAZARDOUS POLYMERIZATION

HAZARDOUS POLYMERIZATION: WILL NOT OCCUR

SECTION 11 - TOXICOLOGICAL INFORMATION

ROUTE OF EXPOSURE

SKIN CONTACT: MAY CAUSE SKIN IRRITATION.

SKIN ABSORPTION: MAY BE HARMFUL IF ABSORBED THROUGH THE SKIN.

EYE CONTACT: MAY CAUSE EYE IRRITATION.

INHALATION: MATERIAL MAY BE IRRITATING TO MUCOUS MEMBRANES AND
UPPER RESPIRATORY TRACT. MAY BE HARMFUL IF INHALED.

INGESTION: MAY BE HARMFUL IF SWALLOWED.

SIGNS AND SYMPTOMS OF EXPOSURE

TO THE BEST OF OUR KNOWLEDGE, THE CHEMICAL, PHYSICAL, AND

TOXICOLOGICAL PROPERTIES HAVE NOT BEEN THOROUGHLY INVESTIGATED.

SECTION 12 - ECOLOGICAL INFORMATION

NO DATA AVAILABLE.

SECTION 13 - DISPOSAL CONSIDERATIONS

APPROPRIATE METHOD OF DISPOSAL OF SUBSTANCE OR PREPARATION

CONTACT A LICENSED PROFESSIONAL WASTE DISPOSAL SERVICE TO DISPOSE
OF THIS MATERIAL. DISSOLVE OR MIX THE MATERIAL WITH A COMBUSTIBLE

SOLVENT AND BURN IN A CHEMICAL INCINERATOR EQUIPPED WITH AN

AFTERBURNER AND SCRUBBER. OBSERVE ALL FEDERAL, STATE, AND LOCAL
ENVIRONMENTAL REGULATIONS.

SECTION 14 - TRANSPORT INFORMATION

DOT

PROPER SHIPPING NAME: NONE

NON-HAZARDOUS FOR TRANSPORT: THIS SUBSTANCE IS
CONSIDERED TO BE NON-HAZARDOUS FOR TRANSPORT.

IATA

NON-HAZARDOUS FOR AIR TRANSPORT: NON-HAZARDOUS FOR AIR
TRANSPORT.

SECTION 15 - REGULATORY INFORMATION

UNITED STATES REGULATORY INFORMATION

SARA LISTED: NO
CANADA REGULATORY INFORMATION
WHMIS CLASSIFICATION: THIS PRODUCT HAS BEEN CLASSIFIED IN
ACCORDANCE WITH THE HAZARD CRITERIA OF THE CPR, AND THE MSDS
CONTAINS ALL THE INFORMATION REQUIRED BY THE CPR.
DSL: YES
NDSL: NO

SECTION 16 - OTHER INFORMATION

DISCLAIMER

FOR R&D USE ONLY. NOT FOR DRUG, HOUSEHOLD OR OTHER USES.

WARRANTY

THE ABOVE INFORMATION IS BELIEVED TO BE CORRECT BUT DOES NOT
PURPORT TO BE ALL INCLUSIVE AND SHALL BE USED ONLY AS A GUIDE. THE
INFORMATION IN THIS DOCUMENT IS BASED ON THE PRESENT STATE OF OUR
KNOWLEDGE AND IS APPLICABLE TO THE PRODUCT WITH REGARD TO
APPROPRIATE SAFETY PRECAUTIONS. IT DOES NOT REPRESENT ANY
GUARANTEE OF THE PROPERTIES OF THE PRODUCT. BIOBASIC,
SHALL NOT BE HELD LIABLE FOR ANY DAMAGE RESULTING FROM HANDLING OR
FROM CONTACT WITH THE ABOVE PRODUCT. SEE REVERSE SIDE OF INVOICE
OR PACKING SLIP FOR ADDITIONAL TERMS AND CONDITIONS OF SALE.
COPYRIGHT 2007 BIOBASIC. LICENSE GRANTED TO MAKE UNLIMITED
PAPER COPIES FOR INTERNAL USE ONLY.