



SAFETY DATA SHEET

[Required under safety and health regulations for shipping and handling]

Version: 2017

Date Updated: September 20, 2017

SECTION 1. ----- PRODUCT AND COMPANY IDENTIFICATION-----

Product Name Ammonium acetate
Product Code(s) ADB0032/AR0032
Recommended Use For Laboratory Research Use Only
 Not for Human or Animal Drug Use

Supplier Bio Basic Inc.
Address 20 Konrad Crescent, Markham, Ontario,
 Canada, L3R 8T4
Telephone (905) 474 4493
Fax (905) 474 5794
For Chemical Emergency Phone# (416) 995 9730

SECTION 2. ----- HAZARDS IDENTIFICATION -----

Emergency Overview

WHMIS Classification

Not Rated

Not a hazardous substance or mixture.

HMIS Classification

Health hazard: 0
Flammability: 0
Physical hazards: 0

Potential Health Effects

Inhalation May be harmful if inhaled. May cause respiratory tract irritation.
Skin May be harmful if absorbed through skin. May cause skin irritation.
Eyes May cause eye irritation.
Ingestion May be harmful if swallowed.

SECTION 3. ----- COMPOSITION/INFORMATION ON INGREDIENTS -----

Chemical Name	EC No.	CAS-No	Weight %
Ammonium acetate	211-162-9	631-61-8	95-100

SECTION 4. ----- FIRST-AID MEASURES-----

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration.

In case of skin contact

Wash off with soap and plenty of water.

In case of eye contact

Flush eyes with water as a precaution.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water.

SECTION 5. ----- FIRE FIGHTING MEASURES -----

Conditions of flammability

Not flammable or combustible.

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special protective equipment for firefighters

Wear self-contained breathing apparatus for firefighting if necessary.

Hazardous combustion products

Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NOx)

Explosion data - sensitivity to mechanical impact

No data available

Explosion data - sensitivity to static discharge

No data available

SECTION 6. ----- ACCIDENTAL RELEASE MEASURES-----

Personal precautions

Avoid dust formation. Avoid breathing vapours, mist or gas.

Environmental precautions

No special environmental precautions required.

Methods and materials for containment and cleaning up

Sweep up and shovel. Keep in suitable, closed containers for disposal.

SECTION 7. ----- HANDLING AND STORAGE-----

Precautions for safe handling

Provide appropriate exhaust ventilation at places where dust is formed.

Conditions for safe storage

Keep container tightly closed in a dry and well-ventilated place.

Recommended storage temperature 2 - 8 °C

Hygroscopic.

SECTION 8. ----- EXPOSURE CONTROLS/PERSONAL PROTECTION-----

Personal protective equipment

Respiratory protection

Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type N95 (US) or type P1 (EN 143) dust masks. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Eye protection

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Skin and body protection

Choose body protection in relation to its type, to the concentration and amount of dangerous substances, and to the specific work-place., The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Hygiene measures

General industrial hygiene practice.

Specific engineering controls

Use mechanical exhaust or laboratory fumehood to avoid exposure.

SECTION 9. - - - - - PHYSICAL AND CHEMICAL PROPERTIES - - - - -**Appearance**

Form	solid
Colour	white

Safety data

pH	6.5 - 7.5 at 77.1 g/l at 25 °C (77 °F)
Melting point/freezing point	Melting point/range: 110 - 112 °C (230 - 234 °F)
Boiling point	Decomposes below the boiling point.
Flash point	No data available
Ignition temperature	No data available
Auto-ignition temperature	No data available
Lower explosion limit	No data available
Upper explosion limit	No data available
Vapour pressure	< 0.001 hPa (< 0.001 mmHg)
Density	1.07 g/cm ³
Water solubility	1,480 g/l at 4 °C (39 °F)
Partition coefficient: n-octanol/water	log Pow: -2.799
Relative vapour density	No data available
Odour	No data available
Odour Threshold	No data available
Evaporation rate	No data available

SECTION 10. - - - - - STABILITY AND REACTIVITY - - - - -**Chemical stability**

Stable under recommended storage conditions.

Possibility of hazardous reactions

No data available

Conditions to avoid

Exposure to moisture

Materials to avoid

Strong oxidizing agents, Strong acids

Hazardous decomposition products

Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NOx) Other decomposition products - No data available

SECTION 11. ----- TOXICOLOGICAL INFORMATION -----

Acute toxicity

Oral LD50

Information given is based on data obtained from similar substances.

Inhalation LC50

No data available

Dermal LD50

Information given is based on data obtained from similar substances.

Other information on acute toxicity

LD50 Intraperitoneal - Mouse - 736 mg/kg

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

Reproductive toxicity

No data available

Teratogenicity

No data available

Specific target organ toxicity - single exposure (Globally Harmonized System)

No data available

Specific target organ toxicity - repeated exposure (Globally Harmonized System)

No data available

Aspiration hazard

No data available

Potential health effects

Inhalation	May be harmful if inhaled. May cause respiratory tract irritation.
Ingestion	May be harmful if swallowed.
Skin	May be harmful if absorbed through skin. May cause skin irritation.
Eyes	May cause eye irritation.

Signs and Symptoms of Exposure

Gastrointestinal disturbance, To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

Synergistic effects

No data available

Additional Information

RTECS: AF3675000

SECTION 12. ----- ECOLOGICAL INFORMATION -----

Toxicity

Toxicity to fish LC50 - Cyprinus carpio (Carp) - 308 mg/l - 48 h
Method: OECD Test Guideline 203

Persistence and degradability

Biodegradability Result: - Readily biodegradable

Bioaccumulative potential

Bioaccumulation is unlikely.

Mobility in soil

No data available

PBT and vPvB assessment

No data available

Other adverse effects

No data available

SECTION 13. ----- DISPOSAL CONSIDERATIONS -----

Product

Offer surplus and non-recyclable solutions to a licensed disposal company.

Contaminated packaging

Dispose of as unused product.

SECTION 14. ----- TRANSPORT INFORMATION -----

IMDG

Not dangerous goods

IATA

Not dangerous goods

SECTION 15. ----- REGULATORY INFORMATION -----

WHMIS Classification

Not Rated

This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations and the MSDS contains all the information required by the Controlled Products Regulations.

SECTION 16. ----- OTHER INFORMATION -----

Issuing Date	09-Feb-2009
Revision Date	20-Sept-2017
Revision Note	No information available.
Recommended Restrictions	No information available

Disclaimer

The information provided on this SDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

End of SDS