

SAFETY DATA SHEET

[Required under safety and health regulations for shipping and handling]

Version: 2021

Date Updated: March 25, 2021

SECTION 1. - - - - - - PRODUCT AND COMPANY IDENTIFICATION - - - - - - -

Product Name ADA buffer Product Code(s) AD0003

Recommended Use For Laboratory Research Use Only

Not for Human or Animal Drug Use

Supplier Bio Basic Inc.

Address 20 Konrad Crescent, Markham, Ontario,

Canada, L3R 8T4

 Telephone
 (905) 474 4493

 Fax
 (905) 474 5794

 For Chemical Emergency Phone#
 (416) 995 9730

SECTION 2. ----- HAZARDS IDENTIFICATION -----

Classification of the substance or mixture

Not a hazardous substance or mixture.

GHS Label elements, including precautionary statements

Not a hazardous substance or mixture

Hazards not otherwise classified (HNOC) or not covered by GHS - none

SECTION 3. - - - - COMPOSITION/INFORMATION ON INGREDIENTS - - - - -

Chemical Name	EC No.	CAS-No	Weight %
N-(2-Acetamido)iminodiacetic acid	247-530-0	26239-55-4	≤100

No components need to be disclosed according to the applicable regulations.

SECTION 4. ----- FIRST-AID MEASURES-----

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration.

In case of skin contact

Wash off with soap and plenty of water.

In case of eye contact

Flush eyes with water as a precaution.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water.

Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (see section 2) and/or in section 11

Indication of any immediate medical attention and special treatment needed

No data available

SECTION 5. ----- FIRE FIGHTING MEASURES -----

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special hazards arising from the substance or mixture

Carbon oxides, Nitrogen oxides (NOx)

Special protective equipment for firefighters

Wear self contained breathing apparatus for fire fighting if necessary.

Hazardous combustion products

Explosion data - sensitivity to mechanical impact

no data available

Explosion data - sensitivity to static discharge

no data available

SECTION 6. ----- ACCIDENTAL RELEASE MEASURES-----

Personal precautions, protective equipment and emergency procedures

Avoid dust formation. Avoid breathing vapours, mist or gas. For personal protection see section 8.

Environmental precautions

Do not let product enter drains.

Methods and materials for containment and cleaning up

Sweep up and shovel. Keep in suitable, closed containers for disposal.

Reference to other sections

For disposal see section 13.

SECTION 7. ----- HANDLING AND STORAGE-----

Precautions for safe handling

Provide appropriate exhaust ventilation at places where dust is formed. Normal measures for preventive fire protection.

For precautions see section 2.

Conditions for safe storage

Keep container tightly closed in a dry and well-ventilated place.

Keep in a dry place.

Storage class (TRGS 510): 13: Non Combustible Solids

Specific end use(s)

Apart from the uses mentioned in section 1 no other specific uses are stipulated

SECTION 8. - - - - EXPOSURE CONTROLS/PERSONAL PROTECTION - - - -

Appropriate engineering controls

General industrial hygiene practice.

Personal protective equipment

Respiratory protection

Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type N95 (US) or type P1 (EN 143) dust masks. Use respirators and components tested and approved

under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Eye protection

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Skin and body protection

Choose body protection in relation to its type, to the concentration and amount of dangerous substances, and to the specific work-place., The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Control of environmental exposure

Do not let product enter drains.

SECTION 9. ----- PHYSICAL AND CHEMICAL PROPERTIES -----

Appearance

Form powder Colour white

Safety data

pH 1.5 - 3.0

Melting point/range: 219 °C (426 °F)

point/freezing point

Boiling point no data available

Flash point no data available Ignition temperature no data available Auto-ignition no data available

temperature

Lower explosion limit no data available
Upper explosion limit no data available
Vapour pressure no data available

Density no data available

Water solubility 9.51 g/l at 20 °C (68 °F)

Partition coefficient: no data available

n-octanol/water

Relative vapour

no data available

density

Odour no data available

Odour Threshold no data available

Evapouration rate no data available

SECTION 10. ------STABILITY AND REACTIVITY -----

Reactivity

No data available

Chemical stability

Stable under recommended storage conditions.

Possibility of hazardous reactions

no data available

Conditions to avoid

No data available

Incompatible Materials

Strong oxidizing agents

Hazardous decomposition products

Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NOx)

Other decomposition products - No data available In

the event of fire: see section 5

SECTION 11. ----- TOXICOLOGICAL INFORMATION -----

Acute toxicity

No data available

Inhalation: No data available Dermal: No data available

No data available

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

no data available

Respiratory or skin sensitisation

no data available

Germ cell mutagenicity

No data available

Carcinogenicity

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as

probable, possible or confirmed human carcinogen by IARC.

Reproductive toxicity

no data available

Teratogenicity

no data available

Specific target organ toxicity - single exposure (Globally Harmonized System)

no data available

Specific target organ toxicity - repeated exposure (Globally Harmonized System)

no data available

no data available
Additional Information RTECS: Not available
SECTION 12 ECOLOGICAL INFORMATION Toxicity No data available
Persistence and degradability no data available
Bioaccumulative potential no data available
Mobility in soil no data available
PBT and vPvB assessment PBT/vPvB assessment not available as chemical safety assessment not required/not conducted
Other adverse effects no data available
SECTION 13 DISPOSAL CONSIDERATIONS
Product Offer surplus and non-recyclable solutions to a licensed disposal company.
Contaminated packaging Dispose of as unused product.
SECTION 14 TRANSPORT INFORMATION
DOT (US) Not dangerous goods
IMDG Not dangerous goods
IATA Not dangerous goods
SECTION 15 REGULATORY INFORMATION
This product has been classified in accordance with the hazard criteria of the Hazardous Products Regulations (HPR) and the SDS contains all the information required by the HPR.
SECTION 16 OTHER INFORMATION
Further information: no limited for paper copy, just for internal uses. For research use only. Not intended for human or animal diagnostic or therapeutic uses.

Aspiration hazard no data available

Synergistic effects

Disclaimer

QF26 Rev 2

5

The information provided on this MSDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

Issuing Date: 25-Mar-2021

End of SDS